




# STEVENS

*Inspired by the care you deliver*



## Medical Supplies and Equipment

# Celebrating Historic Milestones



6 distribution centres in Canada

48 sales representatives across the country to assist you

192 years of medical distribution experience

[www.stevens.ca](http://www.stevens.ca)

 [fb.com/TheStevensCompany](https://fb.com/TheStevensCompany)

 [@StevensCompany](https://twitter.com/StevensCompany)

# Stevens Territory Managers

## (Alternate Care)

My Stevens Account #:



### British Columbia and Yukon

Catrina Herron	Alternate Care—Surrey, Delta, Richmond   LTC—Lower Mainland	catrina.herron@stevens.ca
Cindy Sobkowich	Alternate Care—Langley/White Rock/Tri-cities	cindy.sobkowich@stevens.ca
Clive Gammon	Acute & Alternate Care—North Vancouver Island (VIHA)	clive.gammon@stevens.ca
Dave Jewett	Acute & Alternate Care—South Vancouver Island (VIHA)	dave.jewett@stevens.ca
Grant Lebeau	Acute & Alternate Care—Fraser Valley (FHA)	grant.lebeau@stevens.ca
Scott Hilliard	Acute & Alternate Care—Interior (IHA)	scott.hilliard@stevens.ca
Terry Sayers	Alternate Care—Vancouver (Providence/VCH)	terry.sayers@stevens.ca
Tobin McCallum	Acute & Alternate Care— Lower Mainland/Northern BC/YT (PHSA & NHA)	tobin.mccallum@stevens.ca
Johnny Pashalidis	Acute & Alternate Care— Burnaby/New West/Port Moody (FHA & VCHA)	johnny.pashalidis@stevens.ca



### Midwest

Barbara Olson	Winnipeg Acute Care and LTC	barbara.olson@stevens.ca
Brenda Beavis Frey	Acute & Alternate Care—Northern Saskatchewan	brenda.beavisfrey@stevens.ca
Crystal Sieben	Alternate Care—Edmonton Acute & Alternate Care Northern Alberta	crystal.sieben@stevens.ca
Guy Martin	Acute Care—Edmonton Acute & Alternate Care—Central and Eastern Alberta	guy.martin@stevens.ca
Marcia Brauer	Acute Care—Calgary South Acute & Alternate Care—Southern Alberta	marcia.brauer@stevens.ca
Cindy Frank	Acute, Alternate Care, and LTC—Manitoba	cindy.frank@stevens.ca
Terry Young	Acute & Alternate Care—Southern Saskatchewan	terry.young@stevens.ca
Sarah Maynard	Key Account Manager—Alberta (AHS), NWT	sarah.maynard@stevens.ca



### Ontario

Dylan Chynoweth	Doctors, Clinics, Nursing Homes, Laboratories—Eastern Ontario	dylan.chynoweth@stevens.ca
Don DeSouza	Doctors, Clinics, Nursing Homes, Laboratories— S.W. Ontario (Kitchener/Waterloo to Windsor and Brantford/Hamilton to the Niagara Region)	don.desouza@stevens.ca
Ken Cugno	Hospital Rep/Doctors Clinics/Naturopathic Specialist—Southern Ontario	ken.cugno@stevens.ca
Richard Duguay	Doctors/Clinics/Nursing Homes—Central/Northern Ontario	richard.duguay@stevens.ca
Robert Dimmick	Doctors and Clinics—Mississauga, West Toronto	robert.dimmick@stevens.ca
Terry Schaefer	Doctors and Clinics—Ottawa, Acute Care—Eastern Ontario	terry.schaefer@stevens.ca
Fred McVety	Acute & Alternate Care—Western Ontario	fred.mcvety@stevens.ca
Mark Pagliacci	Acute & Alternate Care—Hamilton/Niagara/Kitchener	mark.pagliacci@stevens.ca
Sylvie Mauro	Alternate Care—Midwives Ontario	sylvie.mauro@stevens.ca



### Québec

Stephanie Lafrance	Alternate Care—Province of Québec	stephanie.lafrance@stevens.ca
Stephanie Picard	Acute & Alternate Care—Eastern Québec	stephanie.picard@stevens.ca
Patrice Marchand	Acute & Alternate Care—Western Québec	patrice.marchand@stevens.ca
Nadia St.Pierre	Acute & Alternate Care—Montreal	nstpierre@stevens.ca



### Atlantic Canada


David MacDonald	Alternate, Acute Care & Veterinary—NB & PEI	david.macdonald@stevens.ca
Daniel Haynes	Alternate & Acute Care—Southern NS	daniel.haynes@stevens.ca
Keith Tobin	Alternate & Acute Care—Northern NS + Veterinary: Nova Scotia	keith.tobin@stevens.ca








# STEVENS

*Inspired by the care you deliver*



## Family Business Since 1830



The Stevens Company has been a pillar in the Canadian healthcare community since 1874 and stands as one of the largest medical supply distributors in Canada today. Originally founded in 1830 in London, England, the company is still owned and operated by the Stevens family, now in their 6<sup>th</sup> generation. Over the decades, Stevens has continued to grow and today has six distribution centres across Canada in Toronto, Vancouver, Calgary, Winnipeg, Halifax, and most recently, Montréal. Stevens has an experienced staff, including 50 sales representatives and product specialists to serve our customers.

**Stevens—Inspired by the care you deliver.**

### British Columbia



19051-27<sup>th</sup> Avenue  
Surrey, BC V3Z 5T1  
Tel: 604-634-3088  
Fax: 604-585-0193  
Toll-free: 1-800-565-8444  
Email: [BCCS@stevens.ca](mailto:BCCS@stevens.ca)

### Manitoba



475 De Baets Street-Unit 2  
Winnipeg, MB R2J 4M3  
Tel: 204-885-9440  
Toll-free fax: 1-888-640-8088  
Toll-free: 1-800-665-0368  
Email: [MWCS@stevens.ca](mailto:MWCS@stevens.ca)

### Atlantic Canada



101 Thornhill Drive-Unit 109  
Dartmouth, NS B3B 1S3  
Tel: 902-468-7582  
Fax: 902-468-7824  
Toll-free: 1-800-565-0765  
Email: [ACCS@stevens.ca](mailto:ACCS@stevens.ca)

### Midwestern Canada



2620-61<sup>st</sup> Avenue S.E.  
Calgary, AB T2C 4V2  
Tel: 403-640-2858  
Fax: 403-640-2976  
Toll-free: 1-800-665-0368  
Email: [MWCS@stevens.ca](mailto:MWCS@stevens.ca)

### Ontario Head Office



425 Rainside Drive  
Brampton, ON L7A 0N8  
Tel: 905-791-8600  
Fax: 905-791-6143  
Toll-free: 1-800-268-0184  
Email: [ONCS@stevens.ca](mailto:ONCS@stevens.ca)

### Quebec



9390 boul des Sciences  
Anjou, QC H1J 3A9  
Tél: 514-352-5652  
Télé.: 514-352-9214  
Sans frais: 1-855-660-7750  
Email: [QCSAC@stevens.ca](mailto:QCSAC@stevens.ca)

# Stevens Repair and Calibration Services



**STEVENS**  
*Inspired by the care you deliver*



Your Sales and Support Company

- After-sales service and yearly calibration
- Manufacturer-trained and certified Repair techs
- Quick turnaround
- Annual Inspections

- Preventative maintenance
- Equipment Audits
- Commissioning
- Training and accreditations

**Estimates  
provided!**



Hillrom™







Equipment / Furniture

7



Diagnostic Equipment and Supplies

16



Surgical / Procedural Instruments

40



Batteries and Bulbs

43



Procedural Supplies

45



Needles and Syringes

52



Infection Control / Sterilization

58



Wound Care

71



IV Supplies

74



Emergency / Resuscitation

77



Paper Products

79




Sundries

83



# STEVENS

*Inspired by the care you deliver*



A full-line medical supplies distributor, The Stevens Company is proud of the companies and products we represent. We are committed to working closely with our valued suppliers to service the ever changing needs of our customer base.

## SOME OF OUR VALUED PARTNERS INCLUDE:



AND MANY, MANY MORE...





Midmark 626 Barrier-Free®  
Examination Chair with the Patient  
Support Rails+ Accessory and  
IQvitals Zone Vital Signs Monitor

## Better BP® is Better Care

Midmark designed the only fully integrated point of care ecosystem to help promote a more consistently accurate blood pressure measurement by targeting 3 key areas using 3 unique products.

### Proper Patient Positioning

The **Midmark® 626 with Patient Support Rails+** is designed to promote AMA-recommended patient positioning for blood pressure capture.<sup>1</sup> The barrier-free low chair height allows most patients to place their feet flat on the floor while the accessory rails support the arm at heart height.

### Accurate, Consistent BP Capture

**IQvitals® Zone™ with the SPRINT BP Protocol** introduces automation at the point of care that can help ensure a higher level of standardization, minimizing human variables while maximizing consistency and data accuracy.

### EMR Connectivity

Seamless connectivity from the **IQvitals Zone** device to the **EMR** saves time and reduces the likelihood of data transcription errors. **Midmark Zone technology** transfers data to the EMR using a secure Bluetooth® Low Energy connection at the point of care.

①  
**Proper Patient  
Positioning**



②  
**Accurate, Consistent  
BP Capture**



③  
**EMR  
Connectivity**

# = Better BP

Learn more at: [midmark.com/stevensBP](https://midmark.com/stevensBP)

# Build Your Clinic with Stevens!



Stevens has everything you'll need to supply your clinic. Whether you're outfitting a new build or simply updating an existing practice, ordering from Stevens is easy.

Simply visit [www.buildmyclinic.ca](http://www.buildmyclinic.ca), browse our product catalogue and add items to your quote. Once you're done, submit your quote and let Stevens do the rest!



[stevens.ca/qr-msec20-1](http://stevens.ca/qr-msec20-1)

Visit Build My Clinic





Ritter

## Ritter® 204 Manual Examination Table

- One-piece seamless drawers
- Stylish and easy to clean
- Removable top system
- Enlarged footstep
- Storage access and volume
- Increased weight capacity

#144-204-011

base only, with drawer

#144-204-012

base only, with drawer heater, and pelvic tilt

#144-002-10138-XX

upholstery seamless top (many colour options available)



stevens.ca/qr-msec20-22

Learn more

Ritter

## Ritter® 224 Barrier-Free® Examination Table

- Gas-spring assisted back
- Power height
- 500lb capacity
- 18"-37" seat height adjustment
- Foot control standard

#144-224-002

base only

#144-002-2009-XX

upholstery top (many colour options available)



stevens.ca/qr-msec20-3

Ritter Barrier-Free Tables

Equipment / Furniture



# Does your exam table, stool, or chair have a tear?

Patch it with **CLEANPatch** or **CLEANPatch-V**.  
The Clinically Validated Solution



stevens.ca/qr-msec20-21

Brewer

## Brewer 5000 Access™ Exam Table

Features the largest, most secure step on the market today that allows patients to safely and easily get on and off the table. The table also offers abundant storage capacity with Brewer's pass-through drawer system and is backed by Brewer's standard 5-year warranty.

- 500lb capacity
- Standard 5-year warranty
- Large, stable four leg step
- Pass-through side drawers eliminates right/left hand setup problems
- Largest storage capacity available on the market today

#120-5000-SP-CN-XX colour choice; see samples on the right



Brewer

## Brewer FLEX Access Exam Table

- Unmatched ergonomic patient access
- Industry leading 700lb weight capacity
- Easy, safe and more independent patient movement
- Multiple options to suit your practice

#120-5700-XX

FLEX Access Power Exam Table with Seamless Top, Pneumatic Backrest

#120-5800-XX

FLEX Access Power Exam Table with Seamless Top, Powered Backrest



stevens.ca/qr-msec20-4

Brewer  
Flex  
Brochure

Brewer

## Brewer 4000 Basic Exam Table

- Simple pneumatic back for easy patient positioning
- Front drawer: 19.8" x 19.2" x 3.8"
- Side drawers: 22.6" x 20.2" x 7"
- 325lb weight capacity
- Standard 3-year warranty
- Easy-to-clean seamless upholstery

#120-4000-XX



## Brewer Upholstery Colour Selection (please specify when ordering)



Clamshell (21)



Gunmetal (22)



Feather (23)



Black Satin (24)



Saddle (25)



Tapestry Red (26)



Cabernet (27)



Cocoa (28)



Deep Sea (29)



Ivy (30)



Deep Sapphire (31)



Blue Fog (32)



Ritter

## Ritter 225 Barrier-Free® Examination Chair

Ritter Barrier-Free examination chairs by Midmark adjust in height from 18-37 inches, improving accessibility to the exam chair and reducing the likelihood of patient or staff injury. Power back.

#144-225-002 base only, with heated drawer, pelvic tilt, and electric receptacle

#144-225-003 base only, with clean assist rollers, heated drawer, and pelvic tilt

#144-002-2009-XX upholstery top (specify colour)



stevens.ca/qr-msec20-5

## Ritter Barrier-Free Tables

Mist-813	Stone-814	Restful Path-858	Iron Ore-851	Deep Earth-854
Lunar Gray-845	Obsidian-857	Oasis-843	Soothing Blue-855	Citrus-849
Cranberry-859	Dark Linen-856	Latte-860	Dream-847	Sandy Retreat-850
Mineral-844	Healing Waters-848	Shaded Garden-853	Robust Brown-852	Curative Copper-846
UltraFree Stone-865	UltraFree Restful Path-868	UltraFree Obsidian-867	UltraFree Harbor-861	UltraFree Latte-870
UltraFree Wheatgrass-862	UltraFree Spice-863	UltraFree Cranberry-869	UltraFree Dark Linen-866	UltraFree Branch-864

Brewer

## Blood Drawing Chair

- Bacteriological, chemical and fungal resistant material
- Easy-to-clean seamless upholstery
- 350lb capacity
- 12 colours available

#120-1500-XX



Ritter

## Ritter® 281 Blood Drawing Chair

- Can enhance your clinical workflow while complementing the style of your space with a range of upholstery colour options
- Quick release mechanism on the arms for simple, efficient height adjustment and inward/outward angle

#144-281-011-XXX



Ritter

## Ritter 280 Side Chair

- 400lb capacity
- 18" seat height
- Contoured, upholstered seat and back section provides comfort

#144-280-001-XX



Ritter

## Ritter 272 Air Lift Stool

- Black composite five-leg base
- 300lb capacity
- Available in standard Ritter upholstery colours (above)
- Without back
- Seat height: 18"-24½"

#144-272-001-XX



Brewer

## One-Step Stool

- Chrome-plated steel
- 350lb capacity
- Non-skid top

#120-11220 with handrail



## Step-on Stool

This safety step-up stool features all-welded steel frame construction, non-tip design, platform, and steel reinforced rubber tips. It also includes a nonskid ribbed rubber mat.

#635-116-460



Brewer

## Square Steel Waste Can

- White enamel finish
- 32 quart/8 gallon capacity
- Dimensions: 11½"D x 12¼"W x 18¼"H

#120-35266



Brewer

## Mayo Instrument Stand

- Stainless steel tray
- Adjustable height
- 2 casters

#120-43460





# Brewer



Nobody protects your practice quite like Brewer.

*With Brewer, you get more than great tables,  
you get the power of our expertise and service.*

**Safeguard patients & staff** with exclusive product innovations that set new standards for safe patient handling.

**Meet more needs in a single bound** with one-stop-shopping convenience for exam, treatment and procedure rooms, including everything from tables, to seating and more.

**Speed turnaround times** with rapid response and pricing support up-front.



**Crush your total cost of ownership** with more value-added features standard and FREE dock-to-dock shipping, including Canada.

**Swing into action quickly** with timely installation help and hands-on customer training.

**Enjoy the ultimate protection** for your power tables that only industry-leading warranties can provide.

*Left to right: Assist Power Procedure Table®.  
Access® High-Low Exam Table.  
FLEX™ Access® High-Low Exam Table.*

*To put the super powers of Brewer and your Safe Patient Handling Specialist to work for you, contact your Stevens Company representative or visit [brewercompany.com](http://brewercompany.com)*



[www.brewercompany.com](http://www.brewercompany.com)

# Health o meter®

P r o f e s s i o n a l

## **500KL-BT**

Eye Level Digital  
Scale with Built-in  
Pelstar Wireless  
Technology

Includes  
new features:  
Pre-Tare, Recall,  
Reweigh and  
Everlock®

Available as  
an accessory  
for new and  
existing scales!



## **500HB**

"Live"


Wrap-Around  
Handlebars for  
500 Series Scales

### **All 500 Series Scales Feature**

- Rust and Dent Proof Polymer Platform
- KG-Only Version Available
- Auto Weight Lock
- 8" Wheelbase for Ease of Movement

Health o meter® Professional also offers pediatric scales, height rods, and large platform scales!

weigh**easier**®



Visit [www.homscales.com](http://www.homscales.com) to learn more about the comprehensive line of Health o meter® Professional Scales.



### Stainless Steel Utility Cart

- 15½" x 24"
  - 300lb capacity
  - 3-shelf
- #152-311



Brewer

### Plastic Utility Cart

- 3 ABS plastic shelves (21¾" x 16")
- 12" between shelves
- 4 non-locking, 3" swivel polyurethane ball-bearing casters
- Total capacity: 200lb
- Easy to assemble
- Dimensions with handle: 25¼"L x 17½"W x 35"H

#120-63530



Brewer

### Two-Hook, Four-Leg IV Pole

- Four-leg, epoxy coated 22" steel base
  - Chrome plated steel tubing
  - Dual wheel nylon casters
  - Easy twist lock height adjustment
  - Height: 51¼"-91¾"
  - 45lb capacity
- #120-11300



Health o meter®

### High-Strength, Wall-Mounted Height Rod

- Height rod measurement range: 55½" (141cm)
  - Measuring tape range: 0"-90" (0cm-230cm)
  - Graduation: +/- 1/16" (1mm)
- #130-205HR



Health o meter®

### Wall Mounted Telescopic Metal Height Rod

- Measuring Range: 2"-90" (6cm-230cm)
  - Graduation: 1/16" (1mm)
  - 2-year limited warranty
- #130-201HR



Health o meter®

### Eye Level Digital Scale

- EMR and Pelstar® wireless connectivity
- Built-in height rod
- 500lb / 220kg capacity
- Body Mass Index (BMI) capability
- 2-year limited warranty

#130-500KL-BT



Health o meter®

### Heavy Duty Eye Level Digital Scale

- EMR connectivity via USB
- Integral digital height rod
- Capacity: 600lb / 272kg
- Resolution: 0.2lb / 0.1kg
- Functions: LB / KG conversion, LB / KG lock out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- 2-year limited warranty

#130-600KL



Health o meter®

### Mechanical Beam Scale

- Capacity: 500lb / 200kg
- Graduation: ¼lb / 100g
- Platform size: 10½"W x 14"D x 3¼"H
- Easy-to-use height rod provides height measurements from 23" to 84"
- Mechanical rotating poise bars allow for easy reading of weight in lb or kg
- 10-year limited warranty

#130-450KL



Choosing the right  
Health-o-Meter  
Scale

stevens.ca/qr-msec20-6

**seca**  
Precision for health

### Seca 216 Mechanical Measuring Rod

- Measuring range: max. 55.3" (138cm) / variable from 1.4"-90" (3.5cm-230cm)
- Graduation length: 1/16" (1mm)
- Weight: 3.7lb (1.7kg)
- Dimensions: 4.7" x 59.1" x 8.5" (120mm x 1,500mm x 217mm)
- Wall mounting

#244-2161814009



**seca**  
Precision for health

### Seca 813 High Capacity Digital Flat Scale

- Capacity: 200kg / 440lb
- Graduation: 100g / 0.2lb
- Power supply: battery
- Measure: 433mm W x 47mm H x 373mm D / 17"W x 1.9"H x 14.7"D
- Net weight: 2.9kg / 6.4lb
- Functions: Auto-HOLD, Step-Off, kg/lbs/sts switch-over, Tap-on automatic switch-on, Automatic switch-off

#244-8131321009



**Health o meter**  
www.healthometer.com

### Mechanical Floor Scale

- 400lb/180kg capacity
- Graduation: 1lb/1kg
- 8" dial (raised 4")
- Platform size: 11"W x 12½"D x 3"H
- 1-year limited warranty

#130-160KL



### Digital Baby Scale

- Capacity: 44lb x 0.5oz (20kg x 10g)
- Convenient lb/kg toggle key
- Large, easy-to-read 1 inch LCD display
- Built-in measuring tape: 0"-22½" (0cm-57cm)
- Recall key to view last weight measurement
- Four AA batteries (not included) or AC adaptor (included)
- Automatic shut-off ensures long battery life
- Two-year limited warranty

#837-107423

**seca**  
Precision for health

### Seca 703 High Capacity Column Scale

- With wireless transmission
- EMR compatible
- Capacity: 300kg / 660lb
- Graduations: 50g < 150kg > 100g / 0.1lb < 330lb > 0.2lb
- Power supply: battery, power adapter
- Measure: 360mm W x 930mm H x 520mm D / 14.2"W x 36.6"H x 20.5"D
- Net weight: 15.2kg / 33.5lb
- Functions: PRINT, Overload protection, Auto-PRINT, Transport casters, TARE, Auto-clear, Pre-TARE, Mother/Child function, HOLD, Auto-HOLD, BMI, kg/lbs/sts switch-over, Automatic switch-off, Wireless interface, SEND, RESET, Acoustic signals can be activated, Damping, Auto-SEND

#244-7031321007



stevens.ca/qr-msec20-7


Seca Scales  
Resource Centre

**seca**  
Precision for health

### Seca® 777 Digital Column Scale w/ Eye Level Display

- Medical-grade, robust scale for years of use with high patient volume
- Built-in measuring rod with 4 - 90" measuring range
- USB interface for EMR integration solutions or PC connectivity
- Large rubber coated castors for excellent movability
- Robust tilt-proof full metal platform, no plastic
- 4-year warranty (double the industry standard)

#244-7771821004



**seca**  
Precision for health

### Seca® 354 Digital Baby Scale

- Baby scale and floor scale for children in one
- Weighing tray easy to detach from base
- Breast-Milk-Intake-Function
- Lightweight and robust for transport
- Capacity: 20kg (44lb)
- Functions: mobile measuring and weighing, TARE, BMIF, auto-HOLD, lbs/kg switch-over, automatic switch-off

#244-3541317004



## General Lights




**Hillrom. WelchAllyn**

### Green Series™ 300 General Exam Light

Perfect combination of intense light output and cool operation, all delivered in a compact design. Ideal for any in-office exam.

#111-44400 with mobile stand  
#111-44410 with table/wall mount

## Procedure Lights



**Hillrom. WelchAllyn**

### Green Series™ 600 Minor Procedure Light

- Ideal for minor procedures and exams in any healthcare environment
- Features 3 LEDs and a wider head specifically for demanding viewing conditions that require exceptionally good light
- Touchless power controls and adjustable light output
- With mobile stand

#111-44600

## Exam Lights

### Heine® EL3 LED Examination Light

- Shadow-free illumination. The compact illumination head (Ø approx. 60 mm) allows convenient, precise positioning, particularly in difficult application situations. Optimal positioning eliminates unwanted shadows in the illuminated field
- Cool-touch closed head design makes disinfection safe and easy
- The light weight head and stable goose neck allow easy and intuitive adjustment and keep the light exactly in the required position
- With wheeled stand



#799-J-008.27.014

### **Hillrom. WelchAllyn** Green Series™ Exam Light IV

- Perfect for the OB/GYN office or for specialties that require very focused light in specific areas
- 3x the typical light output and features a highly flexible light pipe
- Compact, unobtrusive head that ensures pinpoint accuracy and true tissue colour during an exam
- With table/wall mount

#111-48810



### **Ritter** Ritter® 253 LED Exam Light

- ED technology for up to 50,000 hours of useful life
- Designed to consume 45-60% less energy than other lights
- 4,500 K colour temperature
- 1,550 fc intensity (16,500 lux) at 16" focal length
- Sleek, ergonomic and easy-to-clean design

#144-253-016 w/caster stand  
\*table mount options also available







### EDAN M3 Vital Signs Monitor

- Auto/Manual/Continuous/Average BP mode
- Pitch tone
- Real-time measurements display
- Trend table and trend graphs
- Long working-time on battery only
- USB flash disk/Barcode scanner support
- Nurse call function
- LAN/Wi-Fi(optional) connections
- Bi-directional communication with the MFM-CMS central station
- HL7 support via XML files or MFM-CMS central station

#759-302112584

### Hillrom. **WelchAllyn** Connex® Spot Monitor

- SureBP® technology can provide BP in only 15 seconds
- BP averaging to help ensure accurate hypertension diagnosis
- SureTemp® Plus or Braun ThermoScan® PRO 6000 Thermometry
- Choice of leading SpO2 technologies: Masimo®, Nellcor® or Nonin®
- Partnerships with more than 90 leading EMRs including Epic, Cerner and PointClickCare®
- Connex Scoring App calculator for early warning scoring protocols (e.g., NEWS, MEWS)
- Scale connectivity to Detecto, Health-o-meter® and SECA scales

#111-75CE-B SureBP, Covidien SpO<sub>2</sub>, Braun PRO 6000  
 #111-73CT-B SureBP, Nellcor SpO<sub>2</sub>, SureTemp+  
 #111-71WX-B SureBP, Nonin SpO<sub>2</sub>

\*other configurations available



### EDAN M3A Vital Signs Monitor

- Auto/Manual/Continuous/Average BP mode
- Pitch tone
- Real-time measurements display
- Trend table and trend graphs
- Extremely long working-time on battery only
- Barcode scanner support (optional)
- Nurse call function
- LAN/Wi-Fi connections (optional)
- Bi-directional communication with the MFM-CMS central station (optional)

#759-302302012



## The All New Hillrom® Welch Allyn® Spot Vital Signs® 4400 Device

### Simplicity, delivered

You work in a fast-paced environment, so you can't afford to slow down when it comes to capturing the basics. Our next-generation vital signs device is here to help. Experience simple, accurate vitals capture with a solution that fits smoothly into your daily routine. The Spot Vital Signs® 4400 device offers an efficient way to capture, access and document patient vitals so you can focus on what really matters—spending more time with your patients.

### Simple

An intuitive, all-in-one design makes the Spot Vital Signs® device easy to use—for caregivers at any level. Add it to your existing workflow without disruption and begin exams quickly from the touchscreen.

### Accurate

Trusted for clinical accuracy and consistent measurements, the Spot Vital Signs® device can help you make informed treatment decisions at the point of care. The portable design can also help standardize vitals capture throughout your facility and increase productivity for high-volume environments. Plus, blood pressure averaging is now part of the main workflow, enabling you to take and average multiple readings for improved hypertension detection.

### Connected

Save time and help reduce data entry errors by sending patient data to the EMR.

- #111-44XT Spot 4400 with NIBP and SureTemp® Thermometer  
 #111-44WT Spot 4400 with NIBP, SureTemp® Thermometer, and Nonin® SpO2



### Heine ML4 LED HeadLight®

- Up to 90,000 Lux of bright, white light
- Stepless light intensity control
- Binocular loupes and prism optics for video adaption
- Viewing angle adjustable over a wide range
- Multi-coated glass optics for maximum brightness



#799-J-008.31.415

### Bovie® Smoke Shark™ II Smoke Evacuator

- Light-weight, ultra-quiet, and compact design, allowing for simplified operation and storage
- Minimizes exposure to smoke plume during electrosurgical procedures, increasing facility safety and reducing hazardous risk
- The Extended-Life Filter (SF35) offers up to 35 hours of performance, reducing procedure cost by 38%(±) when used at the lowest setting
- 2-year warranty

- #887-SE02 smoke evacuator, without filter or tubing  
 #887-SF35 long-life 35-hr filter  
 #887-SEAS arm stand  
 #887-SETWS-X 7/8" x 6ft tube with wand and tip, sterile



### Heine HR® 2.5x High-Resolution Binocular Loupes

Set includes: optics with i-View on S-Frame, detachable protective lenses, 2 sterilizable swivel levers, retaining cord, cleaning fluid, cleaning cloth, and carrying case. Working distance: 420mm.

#799-C-000.32.356



### Bovie DERM 942™ High Frequency Desiccator

- 40W of power
- Higher resolution in bi-polar mode
- Features that can save you time and money
- Fewer steps needed to accomplish the same procedures

#887-A942



### Bovie® High-Temperature Cauteries

- #887-AA01 fine tip  
 #887-AA03 loop tip  
 #887-AA05 loop tip w/extended 2" shaft

each  
 box of 10  
 each



# Cardiology

Hillrom **WelchAllyn**

## CP 150™ Electrocardiograph

- 7" colour, touch screen display for fast and easy entry of accurate patient information
- Three lead ECG preview
- Optional adult and pediatric MEANS ECG interpretation software offers a second opinion
- Battery operation lets you easily transport your device to the patient
- Store up to 100 test results on the device or transfer to a USB memory stick

#111-CP150-1ENB ECG, non-interpretive, 12-lead



Diagnostic Equipment and Supplies



Hillrom **WelchAllyn**

## 12-Lead, PC-Based Resting ECG

- Wireless ECG connection with the WAM Wireless Acquisition Module.
- Streamlines workflow by capturing and saving an ECG in as few as two clicks.
- Captures best 10 seconds of ECG data from memory, reducing the need for repeat ECGs.
- VERITAS algorithm limits filtering so as not to sacrifice potentially life-saving waveform data.
- Supports DICOM and HL7 communication.
- Application launches from EMR.
- Saves PDFs to any network folder.
- Complies with ACC, AHA and HRS recommendations for adult and pediatric patients.

#111-CC-RXX-AAXX

ECG with AM12 Wired Acquisition Module




## EDAN

### SE-1200 12-Channel ECG

- Large, high resolution colour touch screen
- Alphanumeric keyboard
- One-touch operation
- Enhanced weak signal detection
- Comprehensive filters & anti-noise technology
- Real-time ECG waveforms freezing and review
- Advanced Glasgow algorithm for automatic measurements and interpretation
- Built-in high resolution thermal printer with external USB printer (optional) support
- Extended printing
- Large internal memory which can be enlarged by USB flash disk
- Built-in data management with reports review
- LAN/Wi-Fi (optional) connections
- Bi-directional communication with SE-1515 Data Management System

#759-SE1200EXB







### Aneroid Wall Sphygmomanometer Heine Gamma XXL LF

- Oversized 130mm diameter scale. Can be pivoted for clear, easy reading from any angle
- Latex-free
- 10-year warranty
- Complete with adult cuff

#799-M-000.09.323



### Heine Gamma G7® Aneroid Hand Sphygmomanometer

- Shock proof
- Latex-free
- Exclusive pressure-distribution system
- Pressure-proof manometer
- Micro-filter—protects valve and manometer
- 56mm diameter black-on-white scale is easy to read
- 10-year warranty
- Comes in a zippered pouch

#799-M-000.09.232

with adult cuff

#799-M-000.09.554

family kit (child, adult, and small adult cuffs)



Hillrom. **WelchAllyn**

### Tycos DuraShock Hand Aneroid Sphygmomanometer DS58

- Premium inflation system for comfort
- Unsurpassed reliability in a traditional hand-held design
- Latex-free

#111-DS58-11

with 1 one-tube adult cuff and nylon zipper case

#111-REUSE-09-1TP

one-tube, child cuff (size 9)

#111-REUSE-10-1TP

one-tube, small adult cuff (size 10)

#111-REUSE-11-1TP

one-tube, adult cuff (size 11)

#111-REUSE-12-1TP

one-tube, large adult cuff (size 12)



Hillrom. **WelchAllyn**

### Welch Allyn® 767 Aneroid Sphygmomanometer

- Laser-engraved dial face provides unmatched accuracy
- Recessed dial face and high-contrast red pointer increase visibility and reduce parallax errors
- Lifetime calibration warranty
- Mercury- and latex-free

#111-7670-01

sphyg with reusable adult cuff; wall mount

#111-REUSE-09-2TP

two-tube, child cuff (size 9)

#111-REUSE-10-2TP

two-tube, small adult cuff (size 10)

#111-REUSE-11-2TP

two-tube, adult cuff (size 11)

#111-REUSE-12-2TP

two-tube, large adult cuff (size 12)



Hillrom. **WelchAllyn**

### DS66 Trigger Aneroid Sphygmomanometer, Gold Series

- The DS66 Trigger features an innovative trigger-release design for easy operation
- Shock resistant to maintain accuracy
- Trigger deflation valve for smooth, one-finger operation

This kit includes: DS66 aneroid gauge, durable one-piece adult cuff, and nylon zippered case.

#111-5098-27



### Diagnostix™ 700 Pocket Aneroid Sphyg

- Chrome-plated precision crafted 300mmHg manometer
- Luminescent dial features bold numerals and extended graduation marks for easy reading in all light conditions and from virtually any angle
- Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system prevents mis-cuffing

#334-700-11ABK



Hillrom. **WelchAllyn****Welch Allyn® Green Series™ 777 Integrated Wall System**

- Upgrade to our new Welch Allyn® PanOptic Plus Ophthalmoscope and Welch Allyn® MacroView Plus Otoscope - both featuring bright white LED technology and the option to add your smart device with our iExaminer SmartBracket to save and share exam images when needed.
- Switch your specula to new LumiView™ Clear Single-Use Ear Specula, with up to 8X\* brighter views when compared to a black ear speculum
- Help improve workflow and save energy with handles that automatically turn on/off

#111-777-SM2WAS

Hillrom. **WelchAllyn****Connex® ProBP™ 3400 Digital Blood Pressure Device**

- Small, handheld design makes device easy to carry and capture readings on rounds or in triage settings
- Obtains accurate readings and minimizes error codes, even in the presence of many cases of motion or weak pulse
- USB-wired technology enables connectivity to select EHR/EMR systems, improving documentation accuracy and efficiency
- SureBP™ technology

#111-34XFHT-B

Hillrom. **WelchAllyn****Welch Allyn® 3.5V Diagnostic Set**

Set includes:

- Coaxial LED ophthalmoscope (11710-L)
- MacroView™ Basic LED otoscope (238-2)
- Lithium Ion Rechargeable Handle (71900-USB)
- Hard case

#111-71-SM2LXU

**OMRON****IntelliSense® Digital Blood Pressure Monitor**

- IntelliSense® technology provides automatic cuff inflation, eliminating any need to preset inflation levels or adjust calibration
- Includes four cuff sizes for arm circumference range of 17cm-50cm (7"-20")
- 5-year limited warranty

#786-HEM-907XL monitor

#786-HEM-907-STAND mobile stand with accessory basket

Hillrom. **WelchAllyn****ProBP 2400 Digital Blood Pressure Device**

- 3-Reading Average Mode can give a broader view of patient's true blood pressure
- Irregular heartbeat detection provides data to help you make testing decisions

#111-2400

Hillrom. **WelchAllyn****3.5V Coaxial Ophthalmoscope**

LED illumination produces shadow-free viewing.

#111-11720-L

Hillrom. **WelchAllyn****KleenSpec® Specula**

Universal disposable ear specula for diagnostic otoscopes.

#111-52432-U 2.75mm

#111-52434-U 4.25mm

bag of 850

**Disposable Diagnostic Penlights**

- Sealed concentric beam
- Clearly visible pupil gauge
- Activated by depressing the pocket clip
- Measures 5" long, 1/2" diameter

#001-SC351P

pack of 6



**HEINE** QUALITY  
MADE IN GERMANY



# USE – CLEAN – STAY SAFE

**This is what you can rely on with the HEINE EN 200 wall unit.**

---

When using a HEINE EN 200 Diagnostic Centre, you can count on an easy and safe hygienic reprocessing. Its housing was designed without any openings on the front side and follows the highest safety standards, which stands true for all HEINE products. Specifically, the power supply and the EN 200 itself have several built-in safety protections against electronic malfunctions.

The EN 200 also offers an optimized workflow and excellent energy efficiency. It can be equipped with a number of HEINE instruments, such as otoscopes, ophthalmoscopes and more. Instrument heads available in LED<sup>HO</sup> – eliminating the need for costly and time consuming bulb replacements.





### Heine® EN200 Diagnostic Sets

- Optimized work flow and excellent energy efficiency
- Optimal hygienic reprocessing
- Expandable via USB interface

Contact your local Sales Rep for more configurations.

#799-A-095.12.992B

#799-DIAG-SET-LED-NOBP

#799-DIAG-SET-LED-BP

(Eastern Canada) BETA 200 LED Otoscope and Ophthalmoscope, BP, Ear tip Disp., Wallboard

(Western Canada) BETA 200 LED Otoscope and Ophthalmoscope, Ear tip Disp., Wallboard

BETA 200 LED Otoscope and Ophthalmoscope, BP, Ear tip Disp., Wallboard



### Heine® Mini 3000® Diagnostic Set

Set contains Mini 3000 ophthalmoscope, Mini 3000 F.O. otoscope, 4 reusable tips, AllSpec® disposable tips, and 2 Mini 3000 battery handles with batteries.

#799-D-873.11.021 hard case

#799-D-883.11.021 soft case



### Heine® DELTAone Dermatoscope

- Includes contact plate with scale, USB cord and case
- Field of view of real 21mm (lens 22mm)
- Custom light intensity: 3 degrees of luminous intensity
- 10x magnification
- Toggle function: switch between polarised and non-polarised light at the touch of a button
- No cross-contamination: autoclavable contact plate with scale
- Maintenance-free due to our dustproof design
- Dioptre compensation
- Rechargeable battery with Li-ion technology: no memory effect; you can safely recharge the battery regardless of its charge state
- Wireless connection via iPhone and WLAN with the HEINE Cube System



#799-K-210.28.305

black

#799-K-211.28.305

white

### Heine® DELTA 30 Dermatoscope

- Includes contact plate with scale, USB cord and case
- Largest field of view of real 30 mm (lens 32 mm)
- Crystal-clear image with the new achromatic HEINE optical system
- Excellent colour rendering due to LED HQ
- Digital documentation with the exclusive HEINE DERM app or the HEINE Cube System



#799-K-230.28.305

### AllSpec Oscope Specula

#001-B-000.11.137 4.0mm

#001-B-000.11.138 2.5mm

box of 1000



### Heine® Delta 20 Dermatoscope

- With Beta® 4 USB rechargeable handle
- LED HQ
- Toggle function
- High-resolution optics
- High-performance LEDs

#799-K-262.28.387



Learn more about  
Heine Dermatoscopes



stevens.ca/qr-msec20-8

\*monitor is not included



### PC-Based ECG

- Powerful software with comprehensive functions
- Lead reversal function to easily rectify wrong placement of electrodes
- Auto arrhythmia detection highlights the arrhythmia part on ECG traces automatically
- Event marker enables you to take notes anytime during sampling
- Feature description offers you more confidence for the diagnosis conclusion
- Reports comparison provides you an easy way of comparing ECG traces

#759-SE-1515STR

Hillrom **WelchAllyn**

### ABPM 7100 Blood Pressure Monitor

- Oscillometric technology allows for easy patient hook-up
- Lightweight, compact, and quiet for improved patient compliance
- Dynamic inflation automatically limits maximum cuff pressure for patient comfort
- New sleeve cuffs are machine-washable due to removable bladder
- Latex-free

#111-ABPM-7100S

ABPM 7100 BP monitor with WelchAllyn CardioPerfect® WorkStation software

#111-ABPM-7100 ABPM 7100 BP monitor only



### Tab Electrode ECG Clips

- Red-black
- Latex-free
- PVC-free

#526-EASICLIP

pack of 10



**The Stevens Company has a wide variety of compatible ECG paper.  
Please contact your Stevens Rep for savings today.**

#721-2947455	Schiller AT-2/AT-2 Plus, 210mm x 280mm, red, F-fold sheets	pack of 175	
#721-2847356	Quinton Q3000/4000, no grid, 216mm x 280mm, F-fold sheets	pack of 218	
#721-3584414	Burdick Eclipse+, 216mm x 280mm, red, F-fold sheets	pack of 200	
#721-2197345	PC-Based ECG Graph Paper, red B, 216mm x 280mm sheets	pack of 500	
#721-8422123	Schiller AT-1/SP-1, 90mm x 90mm, red, F-fold sheets	pack of 288	
#721-2384217	Burdick EK-10 & Elite/Elite II, red, 50mm x 45.72m roll	each	
#721-6571572	Marquette Mac/PC, red, 108mm x 23m roll	each	

## ECG Electrodes

Resting			
#800-2360-X	3M Red Dot™ 2cm x 2cm	pack of 100	
#102-31433538-X	Q-Trace™ 5400, popular level of adhesion	pack of 100	
#102-30807732-X	Q-Trace™ Gold 5500, maximum level of adhesion	pack of 100	
#527-1915-100	Fastrace™ 4	pack of 100	
Monitoring			
#800-2560	3M Red Dot™ with foam tape and sticky gel	bag of 50	
#800-2249-50	3M Red Dot™ 6cm, adult, micropore tape backing	bag of 50	
#102-31013926-X	Medi-Trace® 530, non-radiolucent, 1¾" dia. tear drop, foam	pack of 30	
#102-22450	Medi-Trace® 450, non-radiolucent, 1⅜" x 1⅞", foam	case of 1000	
#102-31050522-X	Medi-Trace® 200, non-radiolucent, 1½" dia. tear drop	case of 100	
#527-1700-030	Cleartrace® clear tape, adhesive gel	pack of 30	
#527-1800-030	Suretrace® ECG Electrode	pack of 30	

### Disposable Razor, Single-Blade

#866-75-4008  
box of 100



### 3M™ Red Dot™ Trace Prep

Use to abrade the skin prior to applying electrodes to improve trace quality.

#800-2236

¾" x 196" roll with dispenser







### SureTemp® Plus 690 Electronic Thermometer

- Approximately 10- to 13-second rectal temperature
- Easy-to-use, intuitive operation
- The removable probe well reduces the risk of cross-contamination
- Convenient storage for 25 covers
- A prominent recall button for the display of the last temperature taken
- Ergonomic design with rubber grip

#111-01690-200 thermometer  
#111-05031-105-X disposable probe covers  
**box of 250**



### CareTemp™ Touch-Free Thermometer

- Infrared sensor reads temperature 4-6cm from the center of forehead
- 60 second pulse timer, memory recall button, C/F scale conversion button
- No patient skin contact required, helps minimize the risk of cross-contamination
- Measurement feedback system detects proper operating distance for accuracy
- Large, backlit LCD display that's intuitive and easy to read

#111-105801



### Nellcor™ Portable SpO2 Pulse Oximeter Monitor

- Supports a robust data set, including real-time SpO2 and PR values, SatSeconds alarm management, pleth waveform, blip bar and tabular trend data
- Nellcor™ advanced signal processing technology delivers reliable operation even in difficult monitoring conditions
- LoSat expanded accuracy range (60% to 100% SpO2) when used with Nellcor™ adhesive sensors with OxiMax™ technology
- Compatible with the complete line of Nellcor™ single-patient-use, reusable and specialty sensors with OxiMax™ technology
- 80 hours of data storage capability
- Supports data export to an external personal computer for data analysis and printing functions

#102-PM10N-AD3



### Finger Pulse Oximeter

- Drop-tested durable and moisture-proof
- Perfect spot checks or continuous readings
- Simple one-button operation
- Checks heart rate
- Monitors respiratory function
- Requires two AAA batteries (not included)

#718-7088



### Braun ThermoScan® PRO 6000 Thermometer

- PerfectTemp™ technology provides superior measurement accuracy
- Corrects for shallow placement of probe in ear
- 60 second pulse timer for manual measurement of respiration or pulse
- Integrates with the Connex Spot Monitor to simplify EMR data transfer
- Charging station and recharging battery pack allow charging in just 3-4 hours
- Security tether attaches the thermometer to the cradle

#111-06000-200 thermometer, with small cradles

#111-06000-100 charging station, with battery pack

#111-06000-801 probe covers  
**case of 800**

#111-06000-005-X probe covers  
**box of 200**



### Spectro2™/10 Pulse Oximeter

- Large, bright, easy-to-read LED display for at-a-glance assessment
- Compatible with Docking Station, Attachable Printer, and BCI® spot-check sensor
- Pulse amplitude Index (PI)
- Stores patient data for transfer to PC for improved efficiency and documentation
- Three-year warranty
- Oximeter with adult spot-check finger sensor

#688-WW1000EN



### EDAN H100B Handheld Pulse Oximeter


- Compact & lightweight
- Battery charger stand design (optional)
- LCD display
- Real-time clock
- Pitch tone
- Advanced SpO2 module with equivalent performance as the industry leaders'
- Real-time measurements display
- Plethysmogram display
- Trend table and trend graphs
- Long working-time on battery only
- Large storage capacity
- Data management via PatientCare Viewer software (optional)

#759-H100B





**Note: some models are not available in all of the shown colours**



**3M**

### Littman® Cardiology IV Stethoscope

- Outstanding acoustics with better audibility of high-frequency sounds
- Two sound paths in one tube eliminates rubbing noise of traditional double tubes
- Next-Gen tubing lasts longer and resists stains
- Headset is easily adjusted for individual fit and comfort
- Angled eartubes align with ear canals
- 27" long

#800-6152	black	#800-6155	hunter green
#800-6156	plum	#800-6154	navy blue



Choose the Littman Stethoscope that's right for you



stevens.ca/qr-msec20-9

**3M**

### Littmann® Master Cardiology™ Stethoscope, 27"

A tunable diaphragm and superb craftsmanship offer the ultimate performance in a superior single-sided stethoscope. Meets the exacting standards of healthcare professionals who demand superior acoustic response.



#800-2160	black	#800-2161	black-plated
#800-2163	burgundy	#800-2167	plum
#800-2164	navy blue		

**3M**

### Littmann® Classic II Pediatric and Infant Stethoscopes, 28"

- Solid stainless steel chestpiece, optimally sized for pediatric and infant auscultation
- 1" traditional bell (pediatric) and 3/4" bell (infant) combined with a patented floating diaphragm
- Single-lumen tubing
- Patented 3M® Littmann® Snap Tight soft-sealing eartips
- Comfortably angled, anatomically correct headset
- Patient-friendly nonchill rim and diaphragm



#800-2114	infant, black	#800-2113	pediatric, black
#800-2114R	infant, red	#800-2113R	pediatric, red



### Littmann® CORE Digital Stethoscope

- Connects to Eko software to visualize and share heart sound waveforms
- Up to 40x amplification (at peak frequency, vs. analog mode)
- Active noise cancellation reduces unwanted background sounds
- Toggle between analog and amplified listening modes
- Soft-sealing ear tips provide an excellent acoustic seal and comfortable fit
- Tunable, dual-sided stainless steel chestpiece with open or closed bell

#800-8480-CA



### Littmann® Classic III™ Stethoscope

- Two-sided chestpiece with tunable diaphragms on both the adult and pediatric sides
- Single-piece tunable diaphragm is easy to attach, and easier to clean because its surface is smooth without crevices
- Pediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim
- Next-generation tubing provides longer life due to improved resistance to skin oils and alcohol; less likely to pick up stains
- Small tunable diaphragm is useful for pediatric, small, or thin patients; around bandages; and for carotid assessment
- Stainless steel chestpiece is precision-machined into an aesthetically pleasing, less angular shape. The stem features open side indicator
- Snap-tight eartips have a soft, smooth surface providing a comfortable acoustic seal and comfortable fit
- Five-year warranty
- Length: 69cm / 27.2"

#800-5620 black, standard finish  
 #800-5630 ceil blue, standard finish  
 #800-5832 lavender, standard finish  
 #800-5835 turquoise, standard finish  
 \*more colour options available



### Heine Gamma® C3 Cardio Stethoscope

- Latex-free tubing
- Ideal for cardiology

#799-M-000.09.944



### Adscope® 602 Traditional Cardiology Stethoscope

- Extra deep bell for unsurpassed low frequency response
- Ultra-sensitive diaphragm for greater amplification and crisper high frequency transmission
- 19" tubing
- Lifetime Warranty covers all parts

#334-602BD burgundy  
 #334-602BK black

\*more colour options available



### Adscope® 603 Clinician Stethoscope

- Clinician headset boasts reinforcing yoke molded into flexible PVC tubing and stainless steel binaurals fixed at 15° angle for greater comfort
- 21" tubing
- Lifetime Warranty covers all parts

#334-603BD burgundy  
 #334-603BK black

\*more colour options available











Questions? Email us at [stevens@stevens.ca](mailto:stevens@stevens.ca)




## Choose the right stethoscope to help you perform at your best.

3M<sup>™</sup> Littmann<sup>®</sup> Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the reliable acoustic performance, high quality and comfort you need with a Littmann stethoscope.

Acoustic Performance Levels					
6	7-8	9-10	10+		
					
<b>Lightweight II S.E. Stethoscope</b>	<b>Classic II Pediatric and Infant Stethoscopes</b>	<b>Classic III Monitoring Stethoscope</b>	<b>Cardiology IV Diagnostic Stethoscope</b>	<b>Master Cardiology Stethoscope</b>	<b>CORE Digital Stethoscope</b>
Ideal for taking blood pressure readings and making limited patient assessments.	Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.	Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.	High-performance auscultation. Hear more with amplification, noise cancellation and mobile connectivity.		
Approximate Length	28" / 71 cm	27" / 69 cm	22" / 56 cm, 27" / 69 cm	27" / 68 cm	27" / 69 cm
Approximate Weight	118 g	150 g	167 g, 177 g	185 g	232 g
Limited Warranty	2 years	3 years	7 years	7 years	2 years
Diaphragm	Tunable	Non-Tunable (Floating)	Tunable (2)	Tunable	Tunable (2)
Approximate Chestpiece Diameter	Double-sided 2.0" x 1.8" / 5.2 cm x 4.5 cm 1.4" x 1.3" / 3.5 cm x 3.2 cm	Double-sided Pediatric 1.4" / 3.4 cm 1.1" / 2.7 cm Infant 1.1" / 2.8 cm 0.7" / 1.7 cm	Double-sided Adult side 1.7" / 4.3 cm Pediatric side 1.4" / 3.4 cm	Single-sided 1.9" / 4.8 cm	Double-sided Adult side 1.7" / 4.3 cm Pediatric side 1.4" / 3.4 cm
Features	Teardrop-shaped chestpiece Fits under blood pressure cuff	Specialty designed for the smallest patients	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Up to 40x amplification* Switch from analog listening to digital Connects to Eko software on a smart device** to capture, store and share data

\*At peak frequency, vs. analog mode




The power to learn, in the palm of your hand.





# 3M™ Littmann® Stethoscope Ordering Information and Tube Colour Guide

Components (Chestpiece/Eartube/Stem)											
Lightweight II SE Stethoscope		Classic II Infant Stethoscope	Classic II Pediatric Stethoscope	Classic III™ Monitoring Stethoscope		Cardiology IV™ Diagnostic Stethoscope		Master Cardiology™ Stethoscope	CORE Digital Stethoscope		
Tube Colour											
All Standard Finish	2450 Black	2114 Black	2113 Black	5620 Black	5627 Burgundy	6151 (22") Black	6156 Plum	2160 Black			
	2451 Burgundy	2114R Red	2113R Red	5621 Grey	5633 Pearl Pink	6152 Black	6158 Raspberry	2163 Burgundy			
	2452 Caribbean	2124 Caribbean	2119 Caribbean	5622 Navy	5831 Plum	6154 Navy	6159 Rose Pink	2164 Navy			
	2454 Cell Blue		2122 Raspberry	5623 Caribbean	5832 Lavender	6155 Hunter Green	6184 Burgundy	2167 Plum			
	2456 Pearl Pink			5630 Cell Blue	5835 Turquoise						
All Black				5648 Raspberry	5839 Lemon-Lime						
All Brass Finish				5803 Black	5868 Burgundy	6163 Black	6168 Navy	2161 Black	8480-CA Black		
All Copper Finish								2175 Black			
All Smoke				5809 Chocolate				2176 Black			
				5811 Black							
Rainbow/Brass/Brass		2157 Raspberry	2153 Caribbean	5806 Raspberry	5807 Caribbean						
Brass/Smoke/Smoke											
Rainbow/Black/Black				5870 Black		6165 Black					
Mirror/Smoke/Smoke				5959 Cell Blue		5863 Navy					
Mirror/Standard/Mirror						6177 Black	6170 Burgundy				
Champagne/Smoke/Smoke						6179 Black					
Champagne/Black/Champagne				5861 Black	5864 Burgundy	6176 Burgundy	6190 Caribbean				
High-Polish Rainbow/Black/Black						6242 Navy					
High-Polish Rainbow/Smoke/Smoke						6240 Black	6241 Raspberry				
Coloured Stem Series: (square inset represents stem colour)											
Black/Black/Colour	Standard Finish		Rainbow	6200 Black	6203 Black	6201 Black	6204 Black				
	Black		Mirror	6202 Navy	6232 Black	6238 Grey	6234 Caribbean				
	Brass Finish		Champagne	6239 Plum	6205 Plum						
	Copper Finish		High-Polish Smoke								
	Smoke		High-Polish Rainbow								
High-Polish Smoke/Black/Colour	Standard Finish		Rainbow	5872 Turquoise	5875 Lime						
	Black		Mirror	5873 Grey							
	Brass Finish		Champagne	5960 Plum	5962 Pearl Pink						
	Copper Finish		High-Polish Smoke								
	Smoke		High-Polish Rainbow								



Colours shown are approximate.

1-800-364-3577 | **Littmann.ca** | Also available: Teaching Stethoscope (2138)

© 2021, 3M. All rights reserved. 3M, Littmann, the L Littmann logo, L Littmann Quality and design, Classic III, Cardiology IV, Master Cardiology and the shape of the stethoscope chestpiece are trademarks of 3M. Used under license in Canada. Eko trademarks, including CORE, are marks and/or registered marks of Eko Devices Inc. All other trademarks are property of their respective owners. 2107-21316 E



### ADC® Adscope® 608 Convertible Stethoscope

- Convertible chestpiece provides adult diaphragm/bell or adult/pediatric diaphragm configurations
- Combination diaphragm/bell adult chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability
- ADC's AFD technology provides enhanced acoustic response: light pressure for low, bell-type frequency, firm pressure for high, diaphragm-type frequency
- Extra-large bell (1 $\frac{3}{8}$ " diameter) for unsurpassed low-frequency response
- Ultra-sensitive diaphragm (1 $\frac{3}{4}$ " diameter) for greater amplification and crisper high-frequency response
- Non-chill bell and diaphragm retaining rims for patient comfort
- Clinician headset boasts reinforcing yoke molded into flexible 22" PVC tubing and stainless-steel binaurals fixed at 15° angle for greater comfort

#334-608BK black

\*more colour options available



### ADC® Proscope™ 640 Sprague Stethoscope

- Threaded chestpiece drum of chrome plated zinc permits the use of five interchangeable chestpiece fittings—an adult and pediatric diaphragm, and three bells (1 $\frac{1}{4}$ ", 1",  $\frac{3}{4}$ ")
- Proprietary valve mechanism eliminates acoustic leakage
- Adjustable aluminum binaurals
- Complete with an accessory pouch containing the bell chestpieces, chestpiece key, 2 spare diaphragms, and one pair of silicone eartips
- Traditional 16" double tube configuration

#334-640BK black

\*more colour options available



### ADC® Adscope® 618 Platinum Pediatric Stethoscope

- Distinctive molded-resin animal face chestpiece optimally sized (1 $\frac{1}{2}$ ") for pediatric use
- ADC's AFD technology design switches between low, bell-like frequency and high, diaphragm-like frequency simply by altering pressure
- Non-chill bell and diaphragm retaining rim for patient comfort
- Clinician headset boasts reinforcing yoke molded into the flexible 22" PVC tubing
- Seven interchangeable animal face snap-ons included: Panda, Frog, Koala, Deer, Monkey, Tiger, and Bear

#334-618P pink

#334-618SF seafoam

#334-618W white



### ADC® Adscope® 658 Electronic Stethoscope

- Amplifies sounds up to 18X (24dB)
- Noise-filtering contact microphone eliminates background interference
- Cardiology headset with bi-lumen tubing and stainless-steel binaurals
- Multifunction button controls the frequency mode (bell or diaphragm) and the volume (eight levels)
- Rechargeable Li-ion battery, twenty hours of continuous use on a full charge
- 30" / 175g

#334-658BK black

\*more colour options available

## HUNTLEIGH

### Sonicaid™ FD1 & FD3 Advanced Fetal Doppler

- Advanced real-time auto-correlation fetal heart rate processing
- Active noise reduction for crystal-clear fetal heart sounds
- 2MHz (FD1) & 3MHz (FD3) waterproof transducer suitable for use in a water birthing environment
- Advanced power management with auto switch-off to conserve battery life
- Long robust coiled cable enables optimum transducer position in all circumstances



#684-FD1-P 2MHz  
#684-FD3-P 3MHz

## HUNTLEIGH

### Sonicaid™ D920 Fetal Doppler

- High-sensitivity ultrasound for improved ease of use
- 2MHz waterproof transducer suitable for use in a water birthing environment
- Low battery indicator
- Easy-clean transducer for improved infection control

#684-D920-P


## HUNTLEIGH

### Dopplex® D900 Non-Directional Vascular Doppler

- Clear audio sounds of blood flow or fetal heart
- Enhanced battery management
- Compatible with the complete range of vascular probes
- Large carry bag and headphones included



#684-D900-P doppler  
#684-OP2XS high-sensitivity probe, 2MHz  
#684-OP3XS high-sensitivity probe, 3MHz  
#684-VP5XS high-sensitivity probe, 5MHz  
#684-VP8-HS high-sensitivity probe, 8MHz  
#684-VP10XS high-sensitivity probe, 10MHz



## EDAN

### SonoTrax Fetal Doppler

- A high performance Fetal Doppler with heart rate
- Compact design
- Built-in audio recorder
- Large LCD display & backlight
- Accurate FHR detection with clear sound
- Automatic / Average / Manual fetal heart rate counting
- Highly sensitive, interchangeable, waterproof probes

#759-MS9104660



**hCG Urine Pregnancy Test**

Individually wrapped.

#610-HCG-4.0 **box of 50****One Step hCG Pregnancy Test Cassettes**

Individually wrapped.

#610-ATHCG-U23 **box of 40****Pregnancy Test Cassettes.**

#553-4220-33

**box of 20****BioStrip® hCG Urine Pregnancy Test**

- Simply dip and read
- Highly sensitive, 20mIU/mL
- Bulk package

#553-4220-5 **box of 25****QuickStick® hCG Urine Pregnancy Test**

- A dipstick style test designed for clinical and professional use
- Simply dip and read
- Convenient individual pouch

#553-4220-4K **box of 50****Strep-A® Rapid Tests**

Dipstick-style.

#610-IST-501

**box of 25****Gastrocult® Test Kit**

Two tests in one. Gastrocult is the only test designed specifically for gastric occult blood and pH testing.

#798-66040A

**box of 40****Gastrocult® Applicator Straws**

Two tests in one. Gastrocult is the only test designed specifically for gastric occult blood and pH testing.

#798-66140A

**bag of 40****Gastrocult® Developer**

15mL.

#798-66115A

**QuickVue® Strep A Dipstick Tests**

Kit contains:

- Individually packaged dipsticks (25)
- Sterile throat swabs (25)
- Tubes (25)
- Extraction reagent A: 4M Sodium Nitrite (1)
- Extraction reagent B: 0.2M Acetic Acid (1)
- Positive control: heat inactivated group A streptococcus (1)
- Negative control: heat inactivated group C streptococcus (1)
- Procedure card (1)
- Package insert (1)

#553-4045-20125 **box of 25****Flexible Silver Nitrate Applicators**

- 6", plastic handle
- Silver Nitrate: 75%
- Potassium Nitrate: 25%

#635-118-395 **tube of 100****Hema-Screen™ Occult Blood Tests**

Triple slides. 50 tests.

#726-HSPP-50 **box of 50 x 3****Hemoccult® ICT**

A fecal immunochemical test (FIT) for the detection of bleeding in stool samples.

- No dietary or drug restrictions
- High clinical sensitivity and stability
- 90-day room-temperature storage

#798-395065A

collection cards kit

#798-395067A

test kit


**case of 10****box of 20****Hemoccult® II Dispensapak Patient Screening Kit**

- Most clinically proven
- Trusted for over 40 years
- Most widely used Fecal Occult Blood Test (FOBT)
- Patient education tools and compliance materials

#798-61100A

**case of 300**

#798-61100A-X

**box of 150****Hemoccult® Developer**

15mL

#798-62115







## All strips come in a box of 100

	Bilirubin	Blood in urine	Glucose	Ketone	Leukocytes	Nitrite	pH	Protein	Specific gravity	Urobilinogen
Multistix® 10 SG Urine Test Strips #109-10339493	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Multistix® 8 SG Urine Test Strips #109-10334035		✓	✓	✓	✓	✓	✓	✓	✓	
Multistix® 5 Urine Test Strips #109-10337415		✓	✓		✓	✓		✓		
Uristix® Urine Test Strips #109-10340594			✓					✓		
Chemstrip® 10 Urine Test Strips #549-11203479119	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Chemstrip® 9 Urine Test Strips #549-10398411119	✓	✓	✓	✓	✓	✓	✓	✓		✓
Chemstrip® 7 Urine Test Strips #549-11008552119		✓	✓	✓	✓	✓	✓	✓		
Chemstrip® 5 Urine Test Strips #549-11893467119		✓	✓		✓	✓		✓		
Chemstrip® GP Urine Test Strips #549-10200743119			✓					✓		
BioStrip® U10 Urine Reagent Strips #553-4192-U10	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
BioStrip® U7 Urine Reagent Strips #553-4192-U7		✓	✓	✓	✓	✓	✓	✓		
BioStrip® 2P Urine Reagent Strips #553-4192-2P			✓					✓		

## Drug of Abuse Tests, with Adulterants

### box of 25

	OPI (Opiates)	MDMA (Ecstasy)	MOP (Morphine)	OXY (Oxycodone)	BZO (Oxazepam)	FTY (Fentanyl)	BUP (Buprenorphine)	COC (Cocaine)	AMP (Amphetamine)	EDDP	THC (Cannabis)	HMO
11 Panel Test #610-ATMDU511-1120		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
7 Panel Test #610-ATMDU57-P700			✓	✓	✓	✓		✓	✓	✓		
7 Cup Test #610-ATMDU67-004-C700			✓	✓	✓	✓		✓	✓	✓		
6 Panel Test #610-ATMDU56-P620			✓	✓	✓	✓		✓		✓		
6 Cup Test #610-ATMDU66-004-C635	✓			✓	✓			✓		✓	✓	
5 Panel Test #610-ATMDU55-P500			✓	✓	✓			✓		✓		
5 Cup Test #610-ATMDU65-004-C510			✓	✓	✓			✓		✓		
6 Cup Test #610-ATMDU66-004-C635	✓			✓	✓			✓		✓	✓	

On volume requirements, other configurations available! Talk to your Stevens Representative

### Roche CoaguChek® INRange System

Connected self-testing meter that enables your patients to test their PT/INR whenever, wherever—for confident control of their VKA therapy.

Built-in wireless, Bluetooth® technology allows you to be connected to your patients, who can now transmit their PT/INR results from their meter to you. This will allow you to keep patients on track and in range—anytime, anywhere

#549-08983143001

CoaguChek INRange Kit (does not include strips)

#549-04696522190

CoaguChek control

#549-07671687019

Test Strips **box of 24**



### CoaguChek® INRange system

Connected self-testing meter that enables your patients to test their PT/INR whenever, wherever—for confident control of their VKA therapy.

Built-in wireless, Bluetooth® technology allows you to be connected to your patients, who can now transmit their PT/INR results from their meter to you. This will allow you to keep patients on track and in range—anytime, anywhere.

#549-08983143001

monitor **each**

#549-07671687019

test strips **box of 24**



### Roche Urisys® 1100



- Fast, accurate urine testing with the push of a button
- Chemstrip urine test strips—for visual analysis or for automated testing with the Urisys 1100® urine analyzer
- Unique patented strip design has proven over time to deliver the highest level of performance and accuracy

#549-03617548001

## Our operator-friendly analyzers generate lab-quality results in minutes

### DCA Vantage® Analyzer

A POC system that helps you detect diabetes and prediabetes,\* monitor glycemic control, and detect early kidney disease.

- Simple test procedure requires no sample or reagent preparation.
- Comprehensive onboard tools include patient-trending graphs and a GFR calculator to support effective time-of-visit consultations.
- DCA Systems are used by three out of five physicians who perform HbA1c testing in their office.\*\*

\*GHX Market Intelligence, market data report. 4Q 2014 report.



HbA1c (1 µL blood sample) and ACR (40 µL urine sample) results in 6 and 7 minutes.

### CLINITEK Status+ Urine Chemistry Analyzer

A POC system with a comprehensive urinalysis menu for routine testing, detection of early kidney disease, and hCG pregnancy testing.

- Instrument-read-and-reported urine strip results eliminate the subjectivity of visual interpretation.
- System delivers accurate spot testing on varying urine concentrations with no test preparation.
- Measuring urinary ACR is recommended for screening individuals at risk of kidney disease.†

†2013 NKF/KDIGO Guidelines for Chronic Kidney Disease.



ACR and other essential reagent strip results in 60 seconds.

### Xprecia Stride Coagulation Analyzer<sup>§</sup>

A handheld POC analyzer that helps you monitor patients on Warfarin anticoagulation therapy.

- Interactive touchscreen, step-by-step instructions, and bar-code scanner enhance usability.
- Ergonomic, easy portability, and a small sample volume maximize testing efficiency.
- Dade® Innovin® reagent provides results comparable to laboratory methods.
- Test strip eject button helps minimize potential biohazard exposure.

§Under FDA review. Not available for sale in the U.S. Product availability varies by country.



Results measured across 0.8–8.0 PT/INR range using 6 µL capillary blood sample.

### SIEMENS

#109-10282970 DCA Vantage analyzer

#109-10282131 Optical test cartridge each

#109-10698915 HbA1c reagent

#109-10311480 Microalbumin/Creatine reagent kit box of 10

DCA Vantage full details



stevens.ca/qr-msec20-11

### SIEMENS

#109-10376324 Clinitek Status+ Analyzer

#109-10310482	Chek-Stix control	box of 25
#109-10324060	Azostix	box of 25
#109-10333033	Labstix	box of 100
#109-10340272	Hemastix	box of 50
#109-10340594	Uristix	box of 100
#109-10340700	Albustix	box of 100
#109-10310618	Clinitest hCG	box of 25
#109-10334035	Multistix 8 SG	box of 100
#109-10337415	Multistix 5	box of 100
#109-10339493	Multistix 10 SG	box of 100

### SIEMENS

#109-10714595 Xprecia Stride analyzer kit

#109-10873436 Xprecia PT liquid quality control kit

#109-10714623 Xprecia reagent kit, Stride PT/INR 100 strips (4 x 25)

Xprecia Stride full details



stevens.ca/qr-msec20-10



# Women's Health



## Pro-Medix Disposable Vaginal Specula, Pederson

- For use with light source
- One hand operation
- Adjustable, self-locking
- Colour coded, single use
- Individually wrapped



#001-SC-801S small, blue  
#001-SC-801M medium, green  
#001-SC-801L large, white

**box of 25**



## Pro-Medix Disposable Vaginal Specula, Graves

- One hand operation
- Adjustable, self-locking
- Single use



#001-SC-401S small  
#001-SC-401M medium  
#001-SC-401L large

**bag of 10**



## Welch Allyn KleenSpec® Disposable Vaginal Specula

Premium 590 series.

#111-59000 small  
**box of 24**  
#111-59001 medium  
**box of 24**  
#111-59004 large  
**box of 18**



## Endosampler

- Used for sampling of endometrial tissue
- Unique design allows for superior accuracy
- More reliable sample than traditional endometrial sample devices



#948-022721  
#948-022720 w/10cc syringe



## Welch Allyn KleenSpec® Cordless Illumination System

- Continuous on time of 80 minutes
- Lithium ion battery provides a recharge time of six hours
- Automatic on/off--automatically illuminates in vaginal spectrum, remove it and the light goes off
- Superior light source
- No bulb replacement necessary
- 5,500°K white LED technology for a whiter, brighter light
- Truer tissue rendition and uniform spot—no dark or hot spots—for enhanced visualization of the exam area
- Produces much less heat than traditional lamps for greater comfort

#111-80010



## HealthCare+ Lubricating Jelly

- Water soluble
- Non-conductive
- Bacteriostatic

#248-LUB-035 3.5g pouches  
**box of 145**

#248-LUB-140-X 140g tube  
**each**



## Aquagel® Lubricating Gel

- Bacteriostatic
- Water soluble
- Convenient flip cap
- 142g tube

#457-57-05



## Lubricating Jelly

- Greaseless
- Water soluble
- Bacteriostatic
- Sterile

#002-P903100  
3g, 0.1oz single-use packets  
**box of 144 packets**

#002-P903200 113g/4oz flip-top tube  
**each**





### Vacuum Curettes, Rigid Curved

- High-quality
- Disposable
- Individually packaged

#948-022106 6mm

#948-022112 12mm

\* many other sizes available



### Cytology Brushes

Blue, 7<sup>3</sup>/<sub>4</sub>", plastic handle.

#553-4201-CB-8B-PKG

#553-4201-CB-8B

**pack of 100  
each**



### Cytology Fixative Spray

100mL, non-aerosol.

#553-4240-CYT-017



### OB-GYN/Throat Swabs

- Rayon tip
- 8" paper shaft

#355-57600

**box of 50**



### IUDs

#610-300+ Flexi-T + 300, multi deliveries

#610-380+ Flexi-T + 380, high age risk

#610-300 Flexi-T 300



### Single-Use IUD Insertion & Removal Kit

- Disposable
  - Single use
- Contains:
- Tenaculum Forceps, 10"
  - Littauer Scissors, 8"
  - Uterine Sound, 13"

#190-SM-IUDK-011-312



### Flexible Curettes

- High Quality
- Disposable vacuum curettes
- Flexible
- Dual Port design
- Individually packaged

#948-022006 6mm

#948-022012 12mm

\* many other sizes available

**pack of 50**



### Cervical Scraper

Smoothly polished, high quality, white birch for treating tender tissue. Non-sterile.

#553-4202-51-200-X double end style

**pack of 200**



### Laminaria Cervical Dilator

- 100% natural. Made from a sea-grown plant containing no synthetic materials
- Provides gentle, gradual dilation



Two medium Laminaria before insertion and after 12 hours of insertion.

Before After

- |             |                             |
|-------------|-----------------------------|
| #948-021002 | 2mm, extra small/extra thin |
| #948-021003 | 3mm small                   |
| #948-021004 | 4mm medium                  |
| #948-021005 | 5mm large/thick             |
| #948-021006 | 6mm extra large/extra thick |
| #948-021008 | 8mm jumbo                   |





### Disposable Diagnostic Penlights

- Sealed concentric beam
- Clearly visible pupil gauge
- Activated by depressing the pocket clip
- Measures 5" long, ½" diameter

#001-SC351P

**pack of 6**

### OtoClear®, powered by Aquabot®

Aquabot® is a simple, yet effective pressure spray bottle system that attaches the OtoClear® Ear Irrigation Tips at the end of the wand device, creating a continuous and safe spray for ear irrigation procedures. Healthcare professionals can operate this system with ease through a simple pump and push release motion that offers an uninterrupted spray process.



- |           |   |
|-----------|---|
| #140-7270 | 20 tips, 1 Aquabot® bottle, 1 ear basin |
| #140-7275 | 5 tips, 1 Aquabot® bottle, 1 ear basin  |
| #140-7217 | Aquabot® adaptor wand                   |
| #140-7200 | tips                                    |
| #140-3700 | ear basin                               |

**box of 40  
pack of 3**

✓ LATEX-FREE

### Eclipse® Probe Cover

- Pre-gelled inside with Aquasonic® 100 Ultrasound Transmission Gel
- Non-irritating, non-sensitizing, non-sterile
- Only one probe cover is required for each procedure
- Size: 2.37/1¼"W x 9½"L

#457-38-01

**box of 100**



### Sheathes® Latex Ultrasound Probe Covers

- Rolled
- Individually wrapped
- Non-sterile
- Size: 1¼" dia. x 8" length

#087-10080

**box of 100**



### ThermaSonic® Gel Warmer

- Powerful heating element uniformly heats bottles of gel
- Microprocessor control maintains constant temperature
- "Thermal cutoff" provides protection from overheating
- Low electrical demand
- High-density, impact-resistant composite

- |            |                 |
|------------|-----------------|
| #457-82-01 | holds 1 bottle  |
| #457-82-03 | holds 3 bottles |



### Ultrasound Gel

- Does not contain salt, alcohol, or fragrances
- Water soluble
- Medium viscosity
- Blue

#001-SC-UG250

Stevens Brand, 250mL

#001-SC-UG5000

Stevens Brand, 5L, with 250mL dispenser

#457-01-08

Parker Brand, 250mL

#457-01-50

Parker Brand, 5L, with 250mL dispenser

#457-01-02

Parker Brand, 60g, Doppler® size tube



### Wavelength LV Ultrasound Gel

- Low viscosity

#264-NTPB202X 250mL bottle  
#264-NTPB500X 5L squeezable cube

**case of 24  
each**



### Wavelength CL Ultrasound Gel

- Medium viscosity
- Clear

#264-NTPC502X 5L squeezable bottle, clear w/1 dispenser bottle  
#264-NTPC500X 5L squeezable bottle, clear w/1 dispenser bottle

**case of 4  
each**



### Wavelength TMS Transmission Gel

- Ultra high viscosity

#264-NTPE500X 1L squeezable bottle

**each**



### Wavelength Massage Lotion

#264-NTPM373X 3.8L rigid jug w/one pump

**case of 1**





stevens.ca/qv-msec20-23

Learn more

## Ring Rescue™ Non-Destructive Ring Removal

Ring Rescue™ is simple: apply, inflate, and within seconds the device is working, with uniform and measurable air pressure to decrease finger swelling and therefore decrease finger size, so that the ring can be removed. Return minutes later, remove the device, and attempt to remove the ring intact.

#586-RKM-CA Ring Rescue™ device with lubricant  
#586-RL-CA lubricant refill

## Lighted Ear Curette™

Each box contains 50 tips, 1 light source, and 1 magnification lens.

#140-2210 AngleLoop®, angled 4mm diameter loop  
#140-2220 MicroLoop®, 3mm diameter loop  
#140-2230 WaveCurette®, wave-like 2.5mm wide tip  
#140-2240 FlexLoop®, 4mm wide tip  
#140-2260 InfantScoop®, 2mm wide tip  
#140-2245 variety pack: all 5 tip styles

box



## ALMEDIC Ear Syringe

Automatic 2oz complete ear syringe. Comes with leather case.  
#197-A20-630



## Miltex Vantage Metal Ear Syringe

- Complete with one short bulbous tip, one tapered tip, and shield
- Chrome

#162-V919-381 2oz capacity  
#162-V919-385 4oz capacity



## AMG MEDICAL Aluminum Alloy Tuning Fork

#635-112-110 C-128 (weighted)  
#635-112-114 C-256 (weighted)  
#635-112-118 C-512



## AMG MEDICAL Taylor Percussion Hammer

- Solid rubber head
- Chrome-plated brass handle

#635-112-200



## Miltex Ear Curette Disposable.

#162-19-320 ring tip  
#162-19-321 loop tip  
#162-19-322 small spoon tip  
#162-19-323 large spoon tip  
each



## Elephant Earwasher

- The original model; ideal for removing ear wax
- Cleaning ears is simple and affordable
- Washer and 3 disposable tips

#368-EW earwasher each  
#368-TW spare tips pack of 20



## ALMEDIC Tape Measure

- 60" / 150cm
  - Plastic case with automatic roll-up when not in use
- #197-58-6150



## Microscope Slides

- 3" x 1" (7.6cm x 2.5cm)
  - Single frosted glass
- #024-LAB-034


box of 72



## AMG MEDICAL Queen Square Percussion Hammer

- Flexible nylon handle
- End point for neurology test
- Black rubber ring

#635-112-212 9" handle  
#635-112-218 12" handle





Ask your Stevens Representative about specialized instrument sets from Miltex®



stevens.ca/qr-msec20-13

Detailed product selection by speciality



**Meisterhand**  
(machine-finished)



**Sklar**  
(hand-finished)




**Vantage**  
(floor-grade)

Operating Scissors, 5½",  
Straight, Sharp/Blunt

**S** #444-14-1055

**MH** #162-MH5-16

**V** #162-V95-16



Lister Bandage Scissors, 5½"

**S** #444-11-1055



Multi-Cut Utility Scissors,  
Serrated, 7½"

**S** #444-11-1280



Dressing Forceps, 5", Serrated

**S** #444-19-1050

**MH** #162-MH6-6

**V** #162-V96-6



Adson Dressing Forceps,  
4¾", Serrated

**S** #444-50-3147

**MH** #162-MH6-118

**V** #162-V96-118



Semkin Dressing Forceps, 5", Serrated

**S** #444-22-4150

narrow

**MH** #162-MH6-104

**V** #162-V96-104



Adson Tissue Forceps, 4¾", 1 x 2

**S** #444-50-3047

**MH** #162-MH6-120

**V** #162-V96-120




Halsted Mosquito Forceps, 5", Straight

**S** #444-17-1450

**MH** #162-MH7-2

**V** #162-V97-2



Halsted Mosquito  
Forceps, 5", Curved

**S** #444-17-1550

**MH** #162-MH7-4

**V** #162-V97-4



Splinter Forceps, 4½"

**S** #444-19-3045

fine point

**MH** #162-MH6-304

plain

**V** #162-V96-304

fine serrated



## Schroeder Tenaculum Forceps

<b>S</b>	#444-91-1290	10"
<b>MH</b>	#162-MH30-965	9¾"
<b>V</b>	#162-V930-965	9½"



## Bozeman Dressing Forceps, Double Curved

<b>S</b>	#444-80-1410	10¾"
<b>MH</b>	#162-MH7-622	10½"
<b>V</b>	#162-V97-622	10½"



## Lucae Ear Forceps, 5½", Bayonet, Serrated

<b>S</b>	#444-67-1455
<b>MH</b>	#162-MH19-370
<b>V</b>	#162-V919-370




## Foerster Sponge Forceps, 9½", Straight, Serrated Jaws

<b>S</b>	#444-87-2195
<b>MH</b>	#162-MH7-600
<b>V</b>	#162-V97-600



## Mayo-Hegar Needle Holder

<b>S</b>	#444-20-2052	5¼"
<b>MH</b>	#162-MH8-42	5"
<b>V</b>	#162-V98-42	5"




## Mayo-Hegar Needle Holder, 6"

<b>S</b>	#444-20-2060
<b>MH</b>	#162-MH8-44
<b>V</b>	#162-V98-44




## Baumgartner Needle Holder, 5½", Serrated

<b>S</b>	#444-21-8005
<b>MH</b>	#162-MH8-40
<b>V</b>	#162-V98-40



## Vienna Nasal Speculum, Medium Size

<b>S</b>	#162-20-3
<b>MH</b>	#162-MH20-3
<b>V</b>	#162-V920-3



## Spencer Stitch Scissors, 3½"

<b>S</b>	#444-22-2834
<b>MH</b>	#162-MH9-100
<b>V</b>	#162-V99-100



## Littauer Stitch Scissors


<b>S</b>	#444-23-1190	5½"
<b>MH</b>	#162-MH9-102	4½"
<b>V</b>	#162-V99-102	4½"



## Sterile Disposables


## Iris Scissors, Straight, 4½"

#444-96-2650




## Foerster Sponge Forceps, 9½"

- Straight
  - Serrated
- #444-96-2532




## Adson Forceps

#444-96-2655	serrated, 4¾"
#444-96-2656	1x2 teeth, 4¾"




## Spencer Stitch Scissors, 4½"

#444-96-2596










## Baumgartner Needle Holder, 5½", Serrated.





#444-96-2580










Maxell® Batteries			
#814-723815	AAA size	each	
#814-723443	AA size	each	
#814-721150	9V	each	
#814-723320	C size	pack of 2	
#814-723020	D size	pack of 2	

Heine® Rechargeable Batteries for Beta® Handles			
#799-X-002.99.315	3.5V, NiMH	each	
#799-X-007.99.381	3.5V, Li-ion	each	

Heine® 3.5V Bulbs			
#799-X-002.88.044	for alpha+ Finoff, Glaucotest, F.O laryngoscope handle, Angled F.O. laryngoscope handle, and illumination head for spreadable anal speculum	each	
#799-X-002.88.086	for K180 ophthalmoscope, with AV connector	each	
#799-X-002.88.070	for Beta 200 ophthalmoscope, with AV connector	each	
#799-X-002.88.078	for Beta 200 otoscope, K180 otoscope, multi-purpose illuminator, and Finoff transilluminator	each	

## Welch Allyn® Rechargeable Batteries and Handles

#111-72200	black, 3.5V, NiCd, for 71670 handles	
#111-72000	red, 2.5V, NiCd, for 70700 and 71500 handles	
#111-72600	blue, 2.5V, for PocketScope (72800)	
#111-72300	orange, 3.5V, NiCd, for 71000-A and 71000-C	
#111-71000-A	3.5V, NiCd, direct plug-in rechargeable handle	

Call us for help  
identifying replacement  
parts for your device.

## Welch Allyn® Replacement Bulbs

#111-03000	for 11710, 12400 and 47300 instruments	
#111-04200	for 48700, 48600, 48610, 48625, and 48635 exam lights	
#111-06300	for exam lights (44100 and 44110)	
#111-06500	for MacroView™ (23810 and 23820)	
#111-07800	for 73210, 78000, and 78010 handles	
#111-08800	for 78810 vaginal illuminator	
#111-03800	for PanOptic™ (11800 , 11810, and 11820)	
#111-09500	for Solarc light source (46100 and 49500)	
#111-04700	for standard laryngoscopes	
#111-04900	for ophthalmoscopes (11620, 11630, 11720, 11730, and 11735)	
#111-03100	fits the following instruments: 25020,21700, 20200, 20202, 73500,28100, 27000, 27050, 26538,26530, 41100, and 43300	

## Look™ Sutures

### Silk, Black Braided

#963-781B 4-0, C-6 needle, 18"  
#963-786B 3-0, C-7 needle, 18"

### Nylon


#963-925B 3-0, C-7 needle, 18"  
#963-921B 4-0, C-17 needle, 18"  
#963-922B 4-0, C-6 needle, 18"  
#963-917B 5-0, C-17 needle, 18"  
#963-918B 5-0, C-6 needle, 18"  
#963-916B 6-0, C-17 needle, 18"

### Chromic Gut

#963-560B 3-0, C-6 needle, 27"  
#963-561B 2-0, C-7 needle, 27"  
#963-559B 4-0, C-6 needle, 27"  
#963-1242B 6-0, C-1 needle, 12"

### Plain Gut

#963-555B 2-0, C-7 needle, 27"  
#963-B592N 3-0, C-6 needle, 18", reverse cutting  
#963-554B 3-0, C-6 needle, 27"  
#963-553B 3-0, C-7 needle, 27"  
#963-551B 4-0, C-6 needle, 27"  
**box of 12**



C-1  
11.00mm  
3/8 circle

C-17  
12.00mm  
3/8 circle

C-6  
18.00mm  
3/8 circle

C-7  
24.00mm  
3/8 circle



## Ethicon™ Sutures

### Silk, Black Braided

#652-632G 3-0, X-1 needle, 18"  
#652-683G 4-0, FS-2 needle, 18"  
#652-682G 5-0, FS-2 needle, 18"  
**box of 12**

### Nylon

#652-662G 4-0, FS-2 needle, 18"  
#652-661G 5-0, FS-2 needle, 18"  
#652-663G 3-0, FS-1 needle, 18"  
#652-667G 6-0, C-2 needle, 18"  
**box of 12**


### Plain Gut

#652-612G 3-0, X-1 needle, 18"  
**box of 12**

#652-H810H 3-0, FS-1 needle, 27"  
#652-H822H 3-0, FS-2 needle, 27"  
#652-H821H 4-0, FS-2 needle, 27"  
**box of 36**

### Prolene

#652-8661G 5-0, FS-2 needle, 18"  
#652-8683G 4-0, FS-2 needle, 18"  
#652-8684G 3-0, FS-1 needle, 18"  
#652-8698G 5-0, P-3 needle, 18"  
**box of 12**



C-2  
13.00mm  
3/8 circle

P-3  
13.00mm, multi pass  
3/8 circle

FS-2  
19.00mm  
3/8 circle

X-1  
22.00mm  
1/2 circle

FS-1  
24.00mm  
3/8 circle



### Sofsil™

#102-SS682 5-0, C-13 needle, 18"  
#102-SS632 3-0, C-23 needle, 30"  
**box of 36**

#102-SS684G 3-0, C-14 needle, 18"  
#102-SS683G 4-0, C-13 needle, 18"  
**box of 12**

### Chromic Gut

#102-SG634 5-0, C-13 needle, 30"  
#102-SG635 4-0, C-13 needle, 30"  
#102-SG636 3-0, C-14 needle, 30"  
**box of 36**

### Plain Gut

#102-SG612G 3-0, C-23 needle, 18"  
**box of 12**  
#102-SG810 3-0, C-14 needle, 30"  
#102-SG822 3-0, C-13 needle, 30"  
**box of 36**

### Surgipro™ / Surgipro™ II

#102-SP684G 3-0, C-14 needle, 18"  
#102-SP683G 4-0, C-13 needle, 18"  
#102-SP5698G 5-0, P-13 needle, 18"  
**box of 12**

### Monosof™ Monofilament Nylon

#102-SN663G 3-0, C-14 needle, 18"  
#102-SN662G 4-0, C-13 needle, 18"  
#102-SN661G 5-0, C-13 needle, 18"  
#102-SN667G 6-0, C-1 needle, 18"  
**box of 12**

### Novafil™ Non-Absorbable Monofilament Polybutester

#102-8886442043 3-0, CE-4/C-13 needle, 18"  
#102-8886442033 4-0, CE-4/C-13 needle, 18"  
#102-8886442023 5-0, CE-4/C-13 needle, 18"  
**box of 12**

### Suture Removal Kit

Kit contains:

- 1 x Littauer wire scissors
  - 1 x wire forceps
  - 1 x 4" x 4", gauze
- #626-7100 **kit**



### Suture Removal Kit

Kit contains:

- 1 x Littauer scissors, floor grade
  - 1 x Adson serrated forceps
  - 1 x alcohol prep pad
  - 1 x PVP prep pad
  - 1 x 4" x 4" gauze
- #626-7105





# Many variations of trays and kits available! Ask your Stevens Representative.

## Wrapped Dressing Tray

- Sterile and latex-free.
- Contents (wrapped in a hospital wrap): 1 x 3-compartment tray; 6 x prep balls; 4 x 4" x 4", 8-ply gauze; 2 x plastic forceps; 1 x towel with poly liner; 1 x waste bag with twist tie; and 1 x plastic forceps in the last fold of the hospital wrap.

#626-3042



## Sterile, Latex-Free Laceration Tray, Lidded

Contents:

- 1 x towel
- 1 x drape, fenestrated
- 2 x blue blotting towels
- 1 x Iris scissors
- 1 x Webster needle holder
- 1 x Adson forceps, teeth
- 5 x gauze 8-ply
- 1 x 3-comp tray

#626-7500



✓ LATEX-FREE ✓ STERILE

## Prima® Premium Blade Scalpels

- Disposable
- Stainless-steel
- Comes with safety guard

#560-PSSDS10 #10  
#560-PSSDS11 #11  
#560-PSSDS12 #12  
#560-PSSDS15 #15

box of 10



✓ LATEX-FREE ✓ STERILE

## Miltex® Biopsy Punches

Disposable.

#162-33-31A	1.5mm	#162-33-31	2mm
#162-33-32	3mm	#162-33-33	3.5mm
#162-33-34	4mm	#162-33-35	5mm
#162-33-36	6mm	#162-33-37	8mm



✓ STERILE

## Dynarex Biopsy Punches

The Dynarex Biopsy Punches combine a seamless stainless-steel with a durable, ribbed handle, providing greater comfort and control to the users.

#718-4087	2mm, ivory	#718-4088	2.5mm, green
#718-4089	3mm, grey	#718-4090	3.5mm, yellow
#718-4091	4mm, blue	#718-4093	5mm, pink
#718-4094	6mm, light green	#718-4096	8mm, purple




✓ LATEX-FREE

## Lighted Forceps

Specially designed for foreign body removal. Set includes 1 light source, 1 magnification lens, and box of 10 forceps.

#140-2750



## GluStitch®

GluStitch® is a cyanoacrylate tissue adhesive that is formulated to function as a wound closure tissue adhesive. GluStitch® is indicated for the closure of uncomplicated skin lacerations where there is little or no tissue loss.

#011-201 0.2mL, single-dose vial, clear  
#011-200 0.2mL, single-dose vial, violet



box of 12 each

✓ LATEX-FREE

## Mastisol® Liquid Adhesive

For the most difficult dressing over a protracted period.

#093-3006-02 59ml bottle  
#093-3006-15 15ml bottle



✓ LATEX-FREE

## Detachol® Adhesive Remover

- Non-aerosol
  - Non-explosive
- #093-3007-04 118ml bottle  
#093-3007-15 15ml bottle



## Precise™ Skin Staple Remover

- Scissors-style handle
- Disposable

#800-SR-3 box of 10



## Skin Staple Removal Kit

- 1 x tray
  - 1 x skin staple remover
  - 1 x 3" x 3" gauze sponge
- #261-716-X



✓ LATEX-FREE ✓ STERILE

## Disposable Safety Scalpels

- Stainless steel
- Steadfast disposal lock

#560-SSS-DS-10 #10  
#560-SSS-DS-11 #11  
#560-SSS-DS-15 #15

box of 10



✓ LATEX-FREE ✓ STERILE

## Surgical Scalpel Blades

- Stainless steel

#560-1010 #10  
#560-1011 #11  
#560-1012 #12  
#560-1015 #15

box of 100



✓ LATEX-FREE ✓ STERILE

## Dermal Curettes

Disposable.

#162-33-52	2mm	#162-33-55	5mm
#162-33-53	3mm	#162-33-57	7mm
#162-33-54	4mm		





### Heine UniSpec® Disposable Sigmoidoscopes

Size: 2cm x 25cm.

#799-E-003.18.825 **box of 25**



Hillrom. **WelchAllyn**

### KleenSpec® Disposable Sigmoidoscopes

With obturator. Size: 1.9cm x 25cm.

#111-53130 **box of 25**



### Heine UniSpec® Disposable Anoscopes

Size: 2cm x 8.5cm.

#799-E-003.19.925 **box of 25**



### Sani-Scope® Disposable Anoscopes

Size: 1.8cm x 8.5cm.

#440-82420

**bag of 25**



Hillrom. **WelchAllyn**

### KleenSpec® Disposable Anoscopes

With obturator. Size: 1.9cm x 10cm.

#111-53110

**box of 25**



✓ **STERILE**

### LIDocaine® Local Anesthetic

#389-0122AG01 1% plain, 20mL vial

#389-0122AJ01 1% plain, 50mL vial

#389-0127AG01 2% plain, 20mL vial

#389-0127AJ01 2% plain 50mL vial

#389-0128AG01 2% w/epinephrine, 20mL vial



✓ **STERILE**

### Xylocaine® Jelly, 2%

Lidocaine hydrochloride jelly, USP.

Sterile topical anesthetic.

#114-061 **30mL tube**

✓ **STERILE**

### Xylocaine® Injections

#114-014 1% plain, 50mL

#114-012 1% plain, 20mL

#114-016 1% with epinephrine (1:100M), 20mL

#114-1461 2% with epinephrine (1:200M), 20mL



✓ **STERILE**

### Xylocaine® Lidocaine Hydrochloride Injection USP

▪ Local anesthetic

▪ Sterile

#114-023 2% plain, 50mL

#114-021 2% plain, 20mL





### Utermöhlen® Cryo Professional

- Quickly and effectively remove skin lesions
- Treats various types of skin lesions within 15-40 seconds
- For medical professional use only

#633-613360	set of 30 x 2mm applicators and 30 x 5mm applicators
#633-UTM0169	kit with 1 x 170mL gas canister and 50 x 5mm applicators
#633-UTM0170	kit with 1 x 170mL gas canister and 60 x 2mm applicators
#633-UTM0171	kit with 1 x 170mL gas canister, 30 x 2mm applicators, and 30 x 5mm applicators



### Verruca-Freeze™ Freeze Kit

Non-flammable, odorless alternative to liquid nitrogen that freezes lesions to an effective temperature up to -94°F (-70°C). 3-6 second spray freezes for 20-30 seconds. CryoCones® in 5 convenient sizes and CryoBuds® in 3 sizes limit the freeze to the lesion being treated.

#533-VFK50	complete kit with 150mL canister (50 freezes)	<b>each</b>
#533-VFC65	175mL canister (65 freezes)	<b>each</b>
#533-BUDS	CryoBuds®, small (for 3mm lesions)	<b>pack of 60</b>
#533-BUDM	CryoBuds®, medium (for 6mm lesions)	<b>pack of 60</b>
#533-BUDL	CryoBuds®, large (for 8mm lesions)	<b>pack of 60</b>
#533-TCONES	CryoCones®, soft-tipped reusable 3mm, 5mm, 7mm, 9mm, and 12mm	<b>set of 5</b>



### Histofreezer® Portable Cryosurgical System

- Indicated for 9 types of lesions

#### Kit contains:

- 2 x 80mL canisters
- 24 small (2mm) applicators
- 36 large (5mm) applicators

#933-2879100580



### Miltex® CryoSolutions®

- Rapid freezing (-128°F/-89°C) in a small, self-contained unit
- Autoclavable glass tips
- Controlled temperature helps reduce treatment time
- Cleanable and sterilizable

#162-33510	complete set; includes unit with standard 1mm wide tip, tip protector, one cartridge (23.5g N <sub>2</sub> O), metal pin, user manual, and plastic case
#162-33511	standard, 1mm wide glass tip
#162-33512	2mm wide glass tip
#162-33514	4mm wide glass tip
#162-33515	gynecology tip, non-flexible, 13cm long
#162-33516	10-pack of cartridges (23.5g N <sub>2</sub> O) with threaded safety valve
#162-33517	4-pack of cartridges (23.5g N <sub>2</sub> O) with threaded safety valve



### Powder-Free Stretch Synthetic Vinyl Exam Gloves, Non-Sterile

#001-SC162100 small  
#001-SC162200 medium  
#001-SC162300 large  
#001-SC162350 X-large

**box of 100**



### Powder-Free Synthetic Vinyl Exam Gloves, Non-Sterile

#001-SC161100 small  
#001-SC161200 medium  
#001-SC161300 large  
#001-SC161350 X-large

**box of 100**



### Accelerator-Free Nitrile Exam Gloves, Powder-Free, Non-Sterile, Chemo-Rated

#001-SC363050 X-small  
#001-SC363100 small  
#001-SC363200 medium  
#001-SC363300 large  
#001-SC363350 X-large

**box of 100**




### Nitrile Exam Gloves, Powder-Free, Non-Sterile, Chemo-Rated

#001-SC361050 X-small  
#001-SC361100 small  
#001-SC361200 medium  
#001-SC361300 large  
#001-SC361350 X-large

**box of 200**



Same box.  
Double the  
capacity.



### DermAssist™ Latex, Powder-Free Surgical Gloves

#563-133600 size: 6  
#563-133650 size: 6.5  
#563-133700 size: 7  
#563-133750 size: 7.5  
#563-133800 size: 8  
#563-133850 size: 8.5

**box of 50 pairs**

### Aurelia Vibrant® Latex, Powder-Free Exam Gloves Textured.

862-98225-BX X-small  
862-98226-BX small  
862-98227-BX medium  
862-98228-BX large  
862-98229-BX X-large

**box of 100**



### Pro-Medix Synthetic Vinyl Exam Gloves, Powder-Free

#001-SC160100 small  
#001-SC160200 medium  
#001-SC160300 large  
#001-SC160350 X-large

**box of 150**



### Glove Dispenser

Wire glove box holder,  
epoxy-coated steel.

#001-SC3003 3-box holder  
#001-SC3000 1-box holder





## N95 Masks

AlphaProTech

### Critical Cover® PFL® N-95 Particulate Respirator

- The Magic Arch™ holds the mask away from the nose and mouth for comfortable breathing
- Positive Facial Lock® with the Twist-Seal® chin piece design helps any wearer create a perfect fit with no gaps
- Downward vision is not restricted



#960-695 **box of 35**

3M

### N-95 Particulate Respirator Mask

- Soft inner shell for greater comfort against the face
- May be worn in surgery
- Fluid resistant
- Disposable

#800-1860 regular, teal  
#800-1860S small, teal  
**box of 20**



3M

### Aura™ Health Care N-95 Particulate Respirator and Surgical Mask 1870+

- Fluid resistant and disposable
- Embossed top panel helps reduce eyewear fogging
- White

#800-1870+B Bulk pack  
**case of 440**  
#800-1870+ Box pack  
**box of 20**



## Earloop Masks

✓ LATEX-FREE

### Primagard® 80 Procedure Face Masks

- Ear-loop
- 98% particulate filtration efficiency
- Blue

#263-PG4-1201  
**box of 50**



AlphaProTech

✓ LATEX-FREE

### Critical Cover® AlphaAir® Face Masks

- Outstanding particle filtration down to 0.1 micron
- Full-length aluminum nosepiece helps provide a perfect fit every time
- Blue

#960-BL-5005 **box of 50**



prIME3 ✓ LATEX-FREE

### Primagard® Face Masks

#263-PG4-2001  
Primagard® 80 surgical tie, blue  
**box of 50**



AMD RiteMed

### Distech® Earloop Mask

Blue.  
#355-2115  
**case of 500**



AMD RiteMed

✓ LATEX-FREE

### Medicom® Bouffant Caps, 21"

- Elastic headband for reduced allergic reactions
- Fluid repellant material for optimal protection

#355-8000 white  
#355-8001 blue  
**case of 5 boxes of 100**



### Full Face Shield 10"

- Double-sided Anti-Fog

#001-SC-FFS10  
**case of 200**



### Wooden Tongue Depressors

#001-SC-TDJ-500 junior  
#001-SC-TDA-500 senior  
**box of 500**



### Rayon/Cotton Balls

- Medium
- 0.22g

#589-21176 cotton  
#589-21276 rayon  
**bag of 2000**





# SOL-MILLENNIUM<sup>®</sup>

## Building a Healthier Tomorrow

### SOL-CARE<sup>™</sup> Shield

- Easy one-hand technique of activating protective arm
- Audible click indicates safety mechanism activation



### SOL-CARE<sup>™</sup> Retractable

- Needle retracts into the barrel
- Compatible with all standard hypodermic needles
- Visual and audible confirmation of safety mechanism activation



#### Non-Safety, Sterile, Non-Latex, Single Patient Use, Disposable

#348-180001PP	1mL, 0.01mL GRADUATION, LUER LOCK	box of 100
#348-180011	1mL, 0.01mL GRADUATION SLIP TIP	box of 100
#348-180003	3mL, 0.1mL GRADUATION LUER LOCK	box of 100
#348-180003ST	3mL, 0.1mL GRADUATION SLIP TIP	box of 100
#348-180005	5mL, 0.2mL GRADUATION LUER LOCK	box of 100
#348-180005ST	5mL, 0.2mL GRADUATION SLIP TIP	box of 100
#348-180010	10mL, 0.2mL GRADUATION, LUER LOCK	box of 100
#348-180010ST	10mL, 0.2mL GRADUATION, SLIP TIP	box of 100
#348-180020	20mL, 1mL GRADUATION LUER LOCK	box of 100
#348-180020ST	20mL, 1mL GRADUATION SLIP TIP	box of 100
#348-180030	30mL, 1mL GRADUATION LUER LOCK	box of 50
#348-180060	60mL, 1mL GRADUATION LUER LOCK	box of 30
#348-180060CT	60mL, 1mL GRADUATION CATHETER TIP	box of 30

### SOL-GUARD<sup>™</sup>

- **Transparent Locking Safety Shield with Easy Grip Ridges**—protects during transportation
- **Transport Position**—eliminates the need to recap the needle after draw-up of medication



[stevens.ca/qr-msec20-19](http://stevens.ca/qr-msec20-19)

More sizes and options available!



# SOL-MILLENNIUM<sup>®</sup>

## Building a Healthier Tomorrow

### Sterile, Non-latex, Single-patient Use, Disposable

#### Retractable Safety Syringes Only, Luer Lock

#348-120006IM	SOL-CARE <sup>™</sup> , 3mL
#348-120007IM	SOL-CARE <sup>™</sup> , 5mL
#348-120008IM	SOL-CARE <sup>™</sup> , 10mL
#348-120009IM	SOL-CARE <sup>™</sup> , 20mL (box of 50)

#### Retractable Combos: Safety Syringe With Detachable Needle, Thin Wall Regular Bevel, Luer Lock

#348-100074IM	SOL-CARE <sup>™</sup> , 3mL,	20G x 1",	YELLOW
#348-100084IM	SOL-CARE <sup>™</sup> , 3mL,	20G x 1½",	BLUE
#348-100075IM	SOL-CARE <sup>™</sup> , 3mL,	21G x 1",	GREEN
#348-100072IM	SOL-CARE <sup>™</sup> , 3mL,	21G x 1½",	GREEN
#348-100085IM	SOL-CARE <sup>™</sup> , 3mL,	22G x 1",	BLACK
#348-100076IM	SOL-CARE <sup>™</sup> , 3mL,	22G x 1½",	BLACK
#348-100077IM	SOL-CARE <sup>™</sup> , 3mL,	23G x 1",	BLUE
#348-100078IM	SOL-CARE <sup>™</sup> , 3mL,	25G x ⅝",	ORANGE
#348-100079IM	SOL-CARE <sup>™</sup> , 3mL,	25G x 1",	ORANGE
#348-140072IM	SOL-CARE <sup>™</sup> , 5mL,	21G x 1½",	GREEN
#348-140076IM	SOL-CARE <sup>™</sup> , 5mL,	22G x 1½",	BLACK
#348-160076IM	SOL-CARE <sup>™</sup> , 10mL,	22G x 1½",	BLACK

#### Combos: Syringe With Detachable Safety Needle, Thin Wall Regular Bevel, Luer Lock

#348-12558SN	SOL-CARE <sup>™</sup> Shield, 1mL,	25G x ⅝"	ORANGE, HINGE
#348-12705SN	SOL-CARE <sup>™</sup> Shield, 1mL,	27G x ½"	GRAY HINGED, SHORT BEVEL
#348-13005SN	SOL-CARE <sup>™</sup> Shield, 1mL,	30G x ½"	YELLOW HINGED, SHORT BEVEL
#348-32215SN	SOL-CARE <sup>™</sup> Shield, 3mL,	22G x 1½"	BLACK, HINGED
#348-32310SN	SOL-CARE <sup>™</sup> Shield, 3mL,	23G x 1"	BLUE, HINGED
#348-32510SN	SOL-CARE <sup>™</sup> Shield, 3mL,	25G x 1"	ORANGE, HINGED
#348-32558SN	SOL-CARE <sup>™</sup> Shield, 3mL,	25G x ⅝"	ORANGE, HINGED

#### Safety Needles Only

#348-SN1815	SOL-CARE <sup>™</sup> Shield, 18G x 1½"	PINK THIN WALL, SHORT BEVEL, HINGED
#348-SN2110	SOL-CARE <sup>™</sup> Shield, 21G x 1"	GREEN, THIN WALL REGULAR BEVEL, HINGED
#348-SN2115	SOL-CARE <sup>™</sup> Shield, 21G x 1½"	GREEN, THIN WALL REGULAR BEVEL, HINGED
#348-SN2210	SOL-CARE <sup>™</sup> Shield, 22G x 1"	BLACK, THIN WALL REGULAR BEVEL, HINGED
#348-SN2215	SOL-CARE <sup>™</sup> Shield, 22G x 1½"	BLACK, THIN WALL REGULAR BEVEL, HINGED
#348-SN2310	SOL-CARE <sup>™</sup> Shield, 23G x 1"	BLUE, THIN WALL REGULAR BEVEL, HINGED
#348-SN2558	SOL-CARE <sup>™</sup> Shield, 25G x ⅝"	ORANGE, THIN WALL REGULAR BEVEL, HINGED
#348-SN2510	SOL-CARE <sup>™</sup> Shield, 25G x 1"	ORANGE, THIN WALL REGULAR BEVEL, HINGED
#348-SN2515	SOL-CARE <sup>™</sup> Shield, 25G x 1½"	ORANGE, THIN WALL REGULAR BEVEL, HINGED
#348-SN2712	SOL-CARE <sup>™</sup> Shield, 27G x ½"	GRAY, SHORT BEVEL, HINGED
#348-SN3012	SOL-CARE <sup>™</sup> Shield, 30G x ½"	YELLOW, SHORT BEVEL, HINGED

#### Blunt Fill Needles With and Without Syringe

#348-110022	18G x 1½"	PINK SAFETY, BLUNT FILL
#348-110022F	18G x 1½"	PINK SAFETY WITH 5 MICRON FILTER, BLUNT FILL

#### Needle, Hypodermic, Without Syringe, Thin Wall

#348-111815	18G x 1½",	SHORT BEVEL
#348-112110	21G x 1",	REGULAR BEVEL
#348-112115	21G x 1½",	REGULAR BEVEL
#348-112310	23G x 1",	REGULAR BEVEL
#348-112506	25G x ⅝",	REGULAR BEVEL

#### Needle, Without Syringe, Conventional Wall

#348-111915	19G x 1½",	REGULAR BEVEL
#348-112612	26G x ½",	INTRADERMAL BEVEL
#348-112705	27G x ½",	SHORT BEVEL
#348-113005	30G x ½",	SHORT BEVEL



## Safety Syringes and Needles



### BD® Safety Syringes with Detachable Needles

- Sterile and latex-free
- Safety feature: shielding mechanism
- Luer-Lok™

#333-305789	BD Eclipse™, 1cc, 27G x ½"	<b>box of 50</b>
#333-305780	BD Eclipse™, 1cc, 25G x ⅝"	<b>box of 50</b>
#333-305781	BD Eclipse™, 3cc, 25G x ⅝"	<b>box of 50</b>
#333-305787	BD Eclipse™, 3cc, 25G x 1"	<b>box of 50</b>
#333-305950	allergy treatment tray (25 1cc syringes), BD SafetyGlide™, 27G x ½"	<b>tray</b>

\*other sizes available



### BD® Eclipse™ Hypodermic Needles

- Sterile and latex-free
- Safety feature: shielding mechanism

#333-305766	18G x 1½"
#333-305763	22G x 1½"
#333-305759	25G x ⅝"
#333-305761	25G x 1"
#333-305767	25G x 1½"
#333-305757	30G x ½"

**box of 100**

\*other sizes available



### Terumo® SurGuard®3 Syringes with Safety Needles

- 3 modes of activation—thumb, finger, or hard surface
- Separate finger/thumb grips with ridges and nonslip surfaces
- Hub and needle locks for security post activation
- Broad range of sizes
- Firm injection feel
- Sharpest needles among major brands tested

#691-SG3-01T2516	1cc, 25G x ⅝" needle
#691-SG3-01T2713	1cc, 27G x ½" needle
#691-SG3-03L2238	3cc, 22G x 1½" needle
#691-SG3-03L2325	3cc, 23G x 1" needle
#691-SG3-03L2516	3cc, 25G x ⅝" needle
#691-SG3-03L2525	3cc, 25G x 1" needle

**box of 100**

\*other sizes available



### Terumo® SurGuard®3 Safety Needles

Unique needle locking mechanism to help reduce the risk of needlestick injury.

#691-SG3-2238	22G x 1½" needle
#691-SG3-2516	25G x ⅝" needle
#691-SG3-2525	25G x 1" needle
#691-SG3-2538	25G x 1½" needle
#691-SG3-2713	27G x ½" needle
#691-SG3-3013	30G x ½" needle

**box of 100**

\*other sizes available



## Covidien® Hypodermic Safety Needles

#102-8881850810	Magellan™, 18G x 1"
#102-8881850815	Magellan™, 18G x 1½"
#102-8881850910	Magellan™, 19G x 1"
#102-8881850010	Magellan™, 20G x 1"
#102-8881850015	Magellan™, 20G x 1½"
#102-8881850158	Magellan™, 21G x ⅝"
#102-8881850110	Magellan™, 21G x 1"
#102-8881850115	Magellan™, 21G x 1½"
#102-8881850210	Magellan™, 22G x 1"
#102-8881850215	Magellan™, 22G x 1½"
#102-8881850358	Magellan™, 23G x ⅝"
#102-8881850310	Magellan™, 23G x 1"
#102-8881850558	Magellan™, 25G x ⅝"
#102-8881850510	Magellan™, 25G x 1"

case of 500



# COVIDIEN



## Covidien® Syringe and Safety Needle Combinations

#102-8881511144	0.3cc Monoject™, insulin, 29G x ½"
#102-8881511136	0.5cc Monoject™, insulin, 29G x ½"
#102-8881511110	1cc Monoject™, insulin, 29G x ½"
#102-8881511235	1cc Monoject™, TB, 25G x ⅝"
#102-8881511310	1cc Monoject™, insulin, 30G x 5/16"
#102-8881511336	0.5cc Monoject™, insulin, 30G x 5/16"
#102-8881822215	12cc Magellan™, 22G x 1½"
#102-8881833810	3cc Magellan™, 18G x 1"
#102-8881833558	3cc Magellan™, 25G x ⅝"
#102-8881833510	3cc Magellan™, 25G x 1"

case of 500

case of 500

case of 500

case of 500

case of 500

case of 500

case of 200

case of 400

case of 400

case of 400



## Covidien® Needle Hypo

#102-8881250107	20G x 1.5"
-----------------	------------

case of 1000



## Conventional Syringes and Needles



### BD® Syringe and Needle Combinations

Luer-Lok™ connector.

#333-309623	1cc TB, 27G x ½"	box of 100
#333-309570	3cc, 25G x ⅝"	box of 100
#333-309581	3cc, 25G x 1"	box of 100
#333-309626	1cc TB, 25G x ⅝"	box of 100
#333-305540	allergy treatment tray (25 1cc syringes), regular bevel, 27G x ½"	

\*other sizes available



### BD® PosiFlush™ Saline Syringe

- Designed to eliminate syringe-induced blood reflux
- Greater visibility of syringe contents
- Generates significantly lower pressure compared to standard 3 mL syringes
- Latex-free and preservative-free

#333-306575	S.P., 10mL	box of 30
-------------	------------	-----------



### BD® General Use Syringes

#333-309653	60cc, BD Luer-Lok™	box of 40
#333-302830	20cc, BD Luer-Lok™	box of 48
#333-302995	10cc, BD Luer-Lok™	box of 200
#333-309646	5cc, BD Luer-Lok™	box of 125
#333-309657	3cc, BD Luer-Lok™	box of 200
#333-309628	1cc, BD Luer-Lok™	box of 100
#333-309659	1cc, TB, slip tip	box of 200

\*other sizes available



### BD® Blunt Fill Needles

18G x 1½"

#333-305180

box of 100



### BD® PrecisionGlide® Needles

#333-305195	18G x 1"	#333-305127	25G x 1½"
#333-305196	18G x 1½"	#333-305136	27G x 1¼"
#333-305165	21G x 1"	#333-305109	27G x ½"
#333-305156	22G x 1½"	#333-305106	30G x ½"
#333-305122	25G x ⅝"	#333-305128	30G x 1"
#333-305125	25G x 1"		

box of 100 \*other sizes available



### BD® Ultra-Fine™ II Insulin Syringe & Needle

#333-320440	31G x 5/16"	box of 100
-------------	-------------	------------





### Terumo® Syringe and Needle Combination

#691-SS01T2713 1cc TB, 27G x 1/2"  
 #691-SS03L2525 3cc, 25G x 1"  
 #691-SS03L2516 3cc, 25G x 5/8"


**box of 100** \*other sizes available



### Terumo® Needles

#691-NN2516R 25G x 5/8"  
 #691-NN2525R 25G x 1"  
 #691-NN3013R 30G x 1/2"

**box of 100** \*other sizes available



### Terumo® Conventional Syringes

- Sterile and latex-free
- Single use

#691-SS01T 1cc TB, with luer slip tip  
 #691-SS03L 3cc, luer lock tip  
 #691-SS05L 5cc, luer lock tip  
 #691-SS10L 10cc, luer lock tip  
 #691-SS20L2 20cc, luer lock tip  
 #691-SS30L 30cc, luer lock tip  
 #691-SS60L 60cc, luer lock tip

\*other sizes available

**box of 100**  
**box of 100**  
**box of 100**  
**box of 50**  
**box of 25**  
**box of 25**



### BD® Sharps Collectors

#333-300974 5.1L, counter-balanced door  
 #333-305447 locking wall bracket  
 #333-300440 7.6L, funnel cap, yellow  
 #333-300450 3.1L, inner liner, yellow  
 #333-300460 1.4L, 1-piece, funnel cap, yellow  
 #333-300466 3.1L, 1-piece, funnel cap, yellow  
 #333-367216 1qt, BD Eclipse™ compatible, yellow



### Monoject® Sharps Collectors

#102-8906-X 1 quart, yellow  
 #102-8910 table top holder  
 #102-8506Y-X 5qt, transparent yellow  
 #102-8970Y 2 gallon, transparent yellow, rotor lid  
 #102-8974 table top holder, almond  
 #102-8525Y-X 3 gallon, transparent yellow, hinge lid  
 #102-8529 table top holder, almond



### Daniels EcoShip "Mail Back" Service

**Order it. Use it. Box it. Dial (or email) it. Done!**

- Pickup and delivery to your door
- Assembled sharps containers
- Approved shipping containers
- Certified, approved, compliant, eco-friendly disposal of containers and waste

**Note: sale of this product is restricted in British Columbia.**

#218-D2 2L, pre-assembled **box of 8**  
 #218-D4 4L, pre-assembled **box of 4**



**3M** Science.  
Applied to Life.™

# Put confidence in your hands.

Your partner  
in sterilization  
monitoring solutions.



## Steam sterilization monitoring, now optimized.

Steam  
results in just  
**24**  
minutes\*



Get confidence in just 24 minutes\* – the 3M™ Attest™ 490 Auto-reader gives you fast results and a user-friendly interface. Get the peace of mind that comes from knowing that every patient will receive a high standard of care.

Available in Canada from:

**3M Medical Solutions Division | 3M Canada**  
300 Tartan Drive, London, ON N5V 4M9  
1 800 364-3577 | [www.3M.ca](http://www.3M.ca)

\* When used with 3M™ Attest™ Auto-reader 490 or the 3M™ Attest™ Auto-reader 490H with software version 4.2.7 or greater.

3M, 3M Science, Applied to Life, Attest and Comply are trademarks of 3M.  
Used under license in Canada. Please recycle. © 2020, 3M. All rights reserved.2001-16421-E





Ritter

### M9 and M11 UltraClave® Automatic Sterilizers

- Exclusive automatic-opening door for ease of use and a quick drying time
- Steam-Flush Pressure Pulse Air Removal System ensures proper steam penetration
- Self-program controls enable you to create different cycle parameters to meet specific instrument processing needs
- Pre-programmed controls for convenience

#144-M9-042 9" x 15" chamber  
 #144-M11-042 11" x 18" chamber



Ritter

### M11 & M9 Thermal Printer

This easy-to-install onboard printer documents sterilizer operation for quick and easy record keeping.

#144-9A599001



Ritter

### M11 & M9 Data Logger

Simplify and streamline the documentation of physical and mechanical performance (sterilization time and temperature) of every sterilization cycle.

#144-9A682001



stevens.ca/qr-msec20-14  
 How to operate your M9/M11

Tuttnauer

### EZ Series Automatic Autoclaves w/Printer

- Closed Door Drying maintains sterility and ensures efficient drying of packs and pouches
- High efficiency air pump circulates air in chamber during drying cycle
- Pre-loaded cycles for unwrapped instruments, wrapped/packs, glassware and additional drying cycle.

#594-EZ9P 9" x 18" chamber  
 #594-EZ10P 10" x 19" chamber



midmark

### Midmark® M9D Semi-Automatic Sterilizer

- Manual door
- Self-program controls
- Pre-programmed controls
- LCD display
- Chamber size: 9" x 15"

#144-M9D-042



midmark

### Speed-Clean Sterilizer Cleaner

- Removes hard water build-up and scale deposits from Midmark and Ritter autoclaves and sterilizers
- For use with the M9 UltraClave and M11 UltraClave
- Concentrated

#144-002-0396-05 16oz bottle



3M

### Attest™ 490 Dual Auto-reader

- Intelligent display with time-to-result countdown feature
- Audible alarm sounds for improved monitoring
- Display can be viewed remotely using your smartphone or PC
- Provides clear positive (+) or negative (-) result on LCD display
- For use with 3M™ Attest™ Biological Indicators 1491, 1492V, and 1295\* to provide results in just 24 minutes

#800-490



3M

### Attest® Biological Incubator and Logbook

- Provides reliable incubation of 3M™ Attest™ Biological Indicator for Steam 1261 and 1262
- 15-well capacity
- Dry block pre-set to 56 °C incubation temperature

#800-116 Incubator  
 #800-1266C Logbook



3M

### Attest™ Rapid Readout Biological Indicators

- For 270°F/132°C gravity steam sterilizers
- One-hour readout

#800-1291

box of 50



3M

### Attest® Biological Indicator

Designed for monitoring steam sterilization processes.

#800-1262 box of 100  
 #800-1262P box of 25



3M

### Attest™ Rapid Readout Biological Indicators

- For 250°F/121°C gravity and 270°F/132°C vacuum assisted steam sterilizers
- Three-hour readout

#800-1292

box of 50



3M

### Comply™ SteriGage Steam Chemical Integrator

Class 5 integrating indicator for use in determining whether the sterilization process conditions were met inside each pack. Can be used in all 250-275°F sterilization cycles.

#800-1243B

bag of 100







# STERIKING®

Hospital Packaging Systems


## SEE-THROUGH POUCHES & ROLLS



## ROLL HOLDERS, SEALERS AND ROUTINE CONTROL PRODUCTS FOR SEALING QUALITY



## TRAY LINERS FOR ABSORBING MOISTURE AND FOR EXTRA PROTECTION



More information:  
[www.steriking.com](http://www.steriking.com)

*We care that you pack safely!*

**WIPAK**

## WIPAK RS120A Rotosealer

- Offers practical means for monitoring of sealing parameters
- Temperature is adjustable from 100°C to 200°C, with a tight tolerance of  $\pm 1\%$  only, and its control is in conformance to DIN 58953-7
- The actual temperature is seen on the display

#164-RS120AV



## Premium

#164-R42/31

#164-R41/31

#164-R40/31

#164-R39/31

### WIPAK Flat Sterilization Rolls

#164-R42/31	15cm x 31m (6" x 102ft)
#164-R41/31	10cm x 31m (4" x 102ft)
#164-R40/31	7.5cm x 31m (3" x 102ft)
#164-R39/31	5cm x 31m (2" x 102ft)



### WIPAK Self-Seal Sterilization Pouches

- Standard, with adhesive closure
- Clear
- Flat

#164-SS1	9cm x 20cm (3½" x 8")
#164-SS2	9cm x 25cm (3½" x 10")
#164-SS3	9cm x 57cm (3½" x 22½")
#164-SS4	13cm x 27cm (5" x 10½")
#164-SS5	20cm x 35cm (8" x 14")
#164-SS7	30cm x 45cm (12" x 17½")

**box of 200**

WIPAK Steriking  
Pouches



stevens.ca/qv-msec20-15

## Economy

### Self-Sealing Sterilization Pouches

Featuring dual steam and EtO process indicator arrows that change colour after processing. Pouches are made from medical grade paper and blue tinted transparent film. Triple sealed edges prevent separation or tearing.

#002-P013554	3½" x 5¼"
#002-P013590	3½" x 9"
#002-P015410	5¼" x 10"
#002-P017513	7½" x 13"

**box of 200**



### 3M Comply™ Lead-Free Steam Indicator Tape

An exposure monitor that securely seals packs, and provides a visual assurance that the package has been exposed to the steam sterilization process. The chemical indicator ink, printed in a diagonal striped pattern, changes colour or turns darker when exposed to the sterilization process.

#800-1322-18MM	¾" x 60yd
#800-1322-24MM	1" x 60yd



### Arjo Wiggins® Sterisheets

#855-260-1515	15" x 15"	<b>case of 1000</b>
#855-260-1818	18" x 18"	<b>case of 500</b>
#855-260-2020	20" x 20"	<b>case of 500</b>
#855-260-2424	24" x 24"	<b>case of 500</b>





## Added Assurance For Every Patient.

**Accel® Wipes** deliver effective, one-step cleaning and disinfection with a choice of dwell times.

- Effective against key pathogens – including MRSA, VRE and Norovirus
- Pre-wetted disinfectant cleaner wipes based on proprietary
- Accelerated Hydrogen Peroxide® technology to deliver fast, effective, responsible cleaning performance
- Compatible with most hard, non-porous surfaces
- After use the key ingredient breaks down into oxygen and water

**Accel® INTERvention™**

1 min. dwell time

**Accel® PREvention™**

3-5 min. dwell time



**IntelliCare™ Dispensing System** is a versatile, modern and durable platform, compatible with liquid/gel & foam products.

- Patented built-in technology that allows for touchless & manual activation, and provides low-content alerts
- Robust and reliable ABS-plastic, scratch and UV-resistant
- Dose size switch and optional snap-in placard for ultimate versatility
- Smart packaging that is full product visibility at all times with a patented bottle design that collapses from the back
- Dermatologically tested formulations



**Contact your local Diversey Rep. for more information: 1-800-668-7171**

Accelerated Hydrogen Peroxide® and Design are trademarks of Diversey, Inc.  
All logos and trademarks are owned by or licensed to Diversey, Inc. © 2022 Diversey Inc. All Rights Reserved. 93192-ADV en (22/070)

**[www.diversey.com](http://www.diversey.com)**



Third-party certification  
of environmentally  
preferable products



## Proactive Day-to-Day



## Outbreak and Seasonal



## IntelliCare™ Hand Hygiene Range



### Accelerated Hydrogen Peroxide

- Stable
- Kills germs without being harmful to the user
- Benefits and efficacy have been validated by third party clinical studies
- Designed with cleaning efficiency, germicidal activity, personal health and safety, and the environment in mind



Accel Disinfectant  
Range

stevens.ca/qr-msec20-16



### Accel® PREvention™

Provides the best balance of speed of disinfection, spectrum of disinfection efficacy, cleaning efficacy, health and safety, and environmental preferability.

**EcoLogo certified! 3 Minute General Virucide, Bactericidal, Fungicidal, Tuberculocidal, Broad Spectrum Sanitizer.**

#907-101103731	1L bottle, RTU	each
#907-101103732	4L bottle, RTU	each
#907-100906721	6" x 7" wipes	tub of 160 wipes
#907-100906722	10" x 10" wipes	tub of 60 wipes
#907-100906590-8-X	7½" x 10" wipes, soft-pack	pack of 8



### Accel® INTERvention

Provides superior speed of disinfection, excellent cleaning efficacy, broad spectrum germicidal efficacy, exceptional health and safety profile, environmentally responsible and designed for use in outbreak intervention.

**1 Minute General Virucide, Bactericide, Tuberculocide, 30 Second Sanitizer.**

#907-100906584	10" x 10" wipes	case of 12 x 60
#907-100906585	6" x 7" wipes	tub of 160
#907-101103733	1L bottle, RTU	each
#907-101103736	4L bottle, RTU	each



### IntelliCare™ Dispensing System

- First hybrid dispenser for professional hand care available worldwide, with patented built-in technology that allows for touchless & manual activation, and provides low-content alerts.
- Robust and reliable ABS-plastic, scratch and UV-resistant
- Dose size switch and optional snap-in placard for ultimate versatility

#907-D6205568	Hybrid Dispenser, White
#907-D6205550	Hybrid Dispenser, Black
#907-D6205605	Drip Tray, Gray
#907-D6205592	Drip Tray, Black
#907-100930833	Soft Care® Foam Sanitizer, Pump 536mL
#907-100907877	Soft Care® All Purpose Foam Soap Insert, 1.3L





## Surface Disinfectant Cleaner Max

- Bactericidal, fungicidal, virucidal and tuberculocidal activity
- A 1 minute (2-minute for tuberculosis), ready-to-use disinfectant, decontaminant and cleaner for hospital, medical, surgical and institutional use on inanimate environmental and equipment surfaces

#800-HSD-024-1M 709 mL spray bottle  
 #800-HSW-100-1M 6" x 10" wipes  
 #800-HSW-160-1M 6" x 6¾" wipes  
 #800-HSW-25-1M 6" x 10" wipes  
 #800-HSW-65-1M 9" x 12" wipes



**case of 12**  
**case of 6 tubs of 100**  
**case of 6 tubs of 160**  
**case of 24 packs of 25**  
**tub of 65**

## CaviWipes™ 2.0

- Proven effective against SARS-CoV-2 on hard non-porous surfaces
- Enhanced Efficacy Claims - Effective against over 40 pathogens, including Norovirus and Candida auris for added confidence in preventing healthcare associated infections.
- Excellent Material Compatibility - Compatible with most non-porous, hard surfaces commonly used in clinical settings
- Easy to use label - Makes it quick and easy to find critical information on efficacy, compatibility, and more
- 2-minute universal contact time for bacteria, viruses, and pathogenic fungi (Candida auris)



#517-11-1410-X 6" x 6.75" wipes **160 wipes**

## Pre-Soaks and Cold Sterilization



### Accel® CS 20 Instrument Chemosterilant

Accel CS 20 is designed for use in Healthcare and Personal Services when sterilization in a rapid contact time is required. Accel CS 20 achieves chemical sterilization and high level disinfection in 20 minutes.

- Sporicidal (20 min at 20°C)
- Tuberculocidal (20 min at 20°C)
- Fungicidal (20 min at 20°C)

#909-11405 4L bottle



### EmPower™ Dual Enzyme Instrument Cleaner

Contains two protease enzymes that offer the broadest cleaning action on a variety of protein soils.

#517-10-4100 1 gallon bottle



### MetriClean2™ Cleaner

The alkaline-pH, low-foaming, concentrated cleaner is ideal for stainless steel instruments and cart washers.

#517-10-8100 1 gallon bottle



### Enzyclean® Protease Enzyme Detergent and Pre-Soak

- Neutral pH
- Phosphate-free
- Rust inhibiting
- Low suds

#042-EBL1HC 1 gallon bottle



### Sani-Cloth® Plus Germicidal Disposable Cloth

- Bactericidal
- Tuberculocidal
- Virucidal

#085-Q90172 6" x 6¾" **tub of 160**



### CaviCide® Hardsurface Disinfectant

- Ready to use
- Contains biodegradable detergent
- Bactericidal
- Fungicidal
- Tuberculocidal

#517-11-1000 1 gallon bottle

#517-11-1024 24oz spray bottle



### MetriCide™ 28 High-Level Disinfectant

A long-life, 2.5% buffered glutaraldehyde that can be used and reused as a sterilant and high-level disinfectant for up to 28 days.

#517-11-2801 1 gallon bottle



### BM-28 PLUS™ Activated 2% Glutaraldehyde Disinfectant

Ready for use with no mixing necessary, this shelf-stable solution will disinfect in 10 minutes and sterilize in 10 hours. Contains a corrosion inhibiting agent. Lemon scented. Broad spectrum killing properties. 4L bottle.

#009-28+



### Cidex OPA High-Level Disinfectant

1 gallon bottle.

#421-20394




### Super Germiphene Concentrated Germicidal Detergent and Deodorant

- Bactericidal
- Fungicidal
- Virucidal

#467-SGM-L-CON 454mL bottle





## When Efficacy and Speed Go Hand-in-Hand.

**Give them the power that comes with knowing CaviWipes 2.0 works effectively against viral pathogens of today<sup>2</sup> and tomorrow<sup>3</sup>.**

After 20 years of trusted infection prevention, CaviWipes knows the need for efficacy. With our next-generation CaviWipes 2.0, we're stepping it up with even stronger efficacy claims<sup>4</sup> and a 2-minute contact time. Qualified for the Broad Spectrum Virucidal Claim for all virus types<sup>1</sup>, CaviWipes 2.0 is your single source for infection prevention today<sup>2</sup> and tomorrow<sup>3</sup>. Hands down.



Learn more at  
[www.metrex.com/en-ca/caviwipes-20](http://www.metrex.com/en-ca/caviwipes-20)



**CaviWipes™ 2.0**

1. Health Canada List [www.canada.ca/en/health-canada/services/drugs-health-products/disinfectants/covid-19/list.html](http://www.canada.ca/en/health-canada/services/drugs-health-products/disinfectants/covid-19/list.html) Accessed October 13, 2021

2. On pathogens listed on label.

3. CaviWipes 2.0 has a broad-spectrum virucidal claim. See Health Canada's "Safety and Efficacy Requirements for Surface Disinfectant Drugs 2020" Pg. 9.

4. Compared to CaviWipes™ and CaviWipes1™.



# We tackle the germs, so you can focus on what's at hand.

The PURELL SOLUTION™ is a comprehensive set of products that work together to help you provide better care.

Contact your sales representative to find out how **Canada's #1 brand of hand sanitizer\*** can help keep your practice healthy.





## PURELL® Advanced Hand Rub Foam

- Green certified, fragrance-free, dye-free foam formulation
- Scientifically advanced, 70% formulation outperforms other leading hand sanitizers milliliter for milliliter<sup>1</sup>
- Enhanced with a nourishing blend of four different skin-conditioners that are clinically proven to help maintain skin health<sup>2</sup>
- Compatible with nitrile, vinyl, and latex gloves
- Hypoallergenic

## PURELL® ES4 Push-Style Hand Sanitizer Dispensers

- AT-A-GLANCE™ refills—Monitor product level with just one look
- SANITARY SEALED™ PET refill is easily recycled

#391-5020-01 dispenser **each**  
 #391-5051-02-CAN00 1200mL refill **case of 2**

## PURELL® ES8 Touch-Free Hand Sanitizer Dispensers

- "Energy-on-the-refill" technology—every refill has a built-in energy source, designed to last as long as the refill

#391-7720-01 dispenser **each**  
 #391-7751-02-CAN00 1200mL refill **case of 2**

1. When tested using the ASTM E 1174 Healthcare Personnel Handwash test methodology, against other hand sanitizers commonly used in healthcare facilities, Study # 110103-101, April 5, 2011; # 100907-101, January 6, 2011; 111209-101, March 8, 2012; BioScience Laboratories, Bozeman, MT.  
 2. Kills more germs than other national brands when tested on hands in accordance with the FDA Healthcare Personnel Handwash methodology.



## PURELL® Hand Sanitizing Wipes Alcohol Formula

- Durable, non-linting wipes
- Kills more than 99.99% of most common germs that may cause illness
- Removes light soils and dirt from hands
- Fragrance- and dye-free
- FDA Food Code compliant

#391-9031-06-CAN00  
**case of 6 tubs of 175**



## PURELL® CRT HEALTHY SOAP™

### High Performance Foam

- Remarkably mild formulation removes over 99% of dirt and germs without antibacterial ingredients or harsh preservatives.
- CLEAN RELEASE™ Technology reaches into skin's cracks and crevices 2X better than regular soap to remove dirt and germs<sup>3,4,5</sup>
- A 0.53 log germ removal improvement on dry, irritated skin compared to regular soap—3.4X fewer germs left on skin<sup>6</sup>
- Formulated specifically for dry, sensitive skin and frequently washed hands
- Removes 33% more bio soil than regular soap<sup>7</sup>
- Free of harsh preservatives, antibacterials, parabens, phthalates, fragrance and dye
- Rinses fast and clean for easy gloving—saves an average of 15,000 gallons per year in a 500-bed hospital
- Better wash experience—95% preferred over regular soap<sup>8</sup>

## PURELL® ES4 Push-Style Soap Dispensers

- AT-A-GLANCE™ refills: Monitor product level with just one look
- SANITARY SEALED™ PET refill is easily recycled

#391-5030-01 dispenser **each**  
 #391-5085-02-CAN00 1200mL refill **case of 2**

## PURELL® ES8 Touch-Free Soap Dispensers

- "Energy-on-the-refill" technology—every refill has a built-in energy source, designed to last as long as the refill

#391-7730-01 dispenser **each**  
 #391-7785-02-CAN00 1200mL refill **case of 2**

\* Cleans and moisturizes.

3. Augustine Scientific, Newbury OH, Surface and Interfacial Rheology Analysis, June 30, 2017.  
 4. Bioscience Laboratories, Inc., Bozeman, MT, Study# 170398-101, Evaluation of In-Vivo Germ Removal, July 5, 2017.  
 5. Augustine Scientific, Newbury OH, Ex Vivo Soil Removal Analysis, August 5, 2017.  
 6. GOJO Industries Inc, Akron OH, Evaluation of handwashes for germ removal from ex vivo skin, Study # EV-2017-002, August 15, 2017.  
 7. Augustine Scientific, Newbury OH, Ex Vivo Soil Removal Analysis, August 5, 2017.  
 8. 2017-08-110608 Healthy Soap CRT Foam and Lather Study.





## GOJO® Dispensers

- Easy, professional retrofits
- Ask how your Stevens Rep can assist with installation



True Fit™ Wall Plate

#391-1990-WHT-12 for LTX-12™  
#391-8890-WHT-12 for ADX-12™



### Dispensers and Refills

#391-1920-04-NB	1200mL, LTX-12™ touch-free dispenser, white/white
#391-1903-02-CAN00	PURELL® Advanced Hand Rub, gel, 1200mL refill
#391-8820-06	1200mL, ADX-12™ manual dispenser, white/white
#391-8809-03-CAN00	PURELL® Advanced Moisturizing Foam Hand Rub, 1200mL refill

Messenger™ Dispenser Station  
Requires True Fit™ Wall Plate

#391-1091-WHT-12



Shield™ Floor  
and Wall Protector

#391-1045-WHT-12

## GOJO® Clear & Mild Foam Handwash

- Spa-inspired, green certified foam hand soap that is both fragrance and dye free
- Enriched with moisturizers and skin conditioners
- EcoLogo™ certified
- SMART-FLEX™ bottle design uses 30% less material than HDPE
- Patent-pending controlled collapse technology maintains bottle shape longer while emptying
- Removable pump for easier recycling
- SANITARY SEALED™ refill helps prevent contamination
- Fresh dispensing valve with each refill



### Dispensers and Refills

#391-1984-04-NB	1200mL, LTX-12™ touch-free dispenser, grey/white
#391-1911-02	1200mL refill (foam)
#391-8884-06	1250mL, ADX-12™ manual dispenser, grey/white
#391-8811-03	1250mL refill (foam)
#391-8784-06	700mL, ADX-7™ manual dispenser, grey/white
#391-8711-04	700mL refill (foam)



### PURELL® Advanced Hand Rub

- Kills 99.9% of germs on hands, including *Staphylococcus aureus*, *Escherichia coli* and *Staphylococcus epidermis*<sup>1</sup>
- Fragrance-free, dye-free gel formulation
- Scientifically advanced, 70% formulation outperforms other leading hand sanitizers milliliter for milliliter<sup>2</sup>



#391-9651-24-CAN00-X	118mL bottle, flip-top
#391-9652-12-CAN00-X	236mL pump bottle

1. NPN 80047004. PURELL Advanced Hand Rub gel formulation.  
2. When tested using the ASTM E 1174 Healthcare Personnel Handwash test methodology, against other hand sanitizers commonly used in healthcare facilities, Study # 110103-101, April 5, 2011; # 100907-101, January 6, 2011; 111209-101, March 8, 2012; BioScience Laboratories, Bozeman, MT.

## GOJO® Infection Prevention Video



stevens.ca/qr-msec20-17



### Isagel® Hand Sanitizing Gel

- Added moisturizers
  - CHG and latex glove compatible
- #310-1644 115mL bottle  
#310-1645 621mL bottle with pump



# Instrument Cleaning Guide

## Spaulding's Classification of Medical Equipment/Devices and Required Level of Processing/Reprocessing

Classification	Definition	Level of Processing/Reprocessing
Critical Equipment/Device	Equipment/Device that enters sterile tissues, including the vascular system	Cleaning followed by Sterilization
Semi-critical Equipment/Device	Equipment/Device that comes in contact with non-intact skin or mucous membranes but does not penetrate them	Cleaning followed by High Level Disinfectant (as a minimum requirement). Sterilization is preferred
Noncritical Equipment/Device	Equipment/Device that touches only intact skin and not mucous membranes, or does not directly touch the client/patient/resident	Cleaning followed by Low Level Disinfectant (in some cases, cleaning alone is acceptable)

## Cleaning



Never keep used/dirty instruments in a dry container. Make sure BLOOD and DEBRIS don't have a chance to dry. If you are unable to process them right away use:

An Enzymatic Cleaner to help to remove organic material  
 #042-EBL1HC Enzyclean® instrument cleaner, pre-soak single, 1 gallon  
 #517-10-4100 EmPower™ dual enzyme instrument cleaner, 1 gallon

Or a neutral pH detergent:  
 #042-B1HC Low suds instrument cleaner, 1 gallon

**Always rinse off chemicals with distilled water.**



**ULTRASONIC CLEANING** is usually recommended as the best way to clean surgical instruments.  
 #594-CSU1 Ultrasonic cleaner, 1 gallon

**MANUAL CLEANING:** Use stiff brushes to make sure all visible dirt is removed.  
 #444-10-1651-X Nylon bristles instrument brush



Use a surgical instrument cleaner  
 #071-002900 Organisol instrument cleaner, 4lb container, with scoop

**Make sure everything looks sharp and all moving parts are working correctly.**



## High Level Disinfection



Keep in mind safety, efficacy, compatibility, and the environment.

### Accelerated Hydrogen Peroxide

Not recommended for instruments with brass, copper, or tungsten carbide.

**Prevention—5-minute HLD for respiratory and plastics. 6-hour chemosterilant**  
 #909-11505 4L bottle

**CS20—20-minute chemosterilant**  
 #909-11405 4L bottle



## Sterilization



It is recommended that all surgical instruments be autoclaved, if possible.  
#144-M11-042 Midmark M11 autoclave, 11" x 18" chamber

Always use distilled water in your autoclave.

#781-AQ00334 Distilled water, 4L, case of 4 (Western Canada Only)

#320-72404 Distilled water, 4L, case of 4 (Ontario, Quebec and Eastern Canada Only)



Lubricate all moving parts.

#162-3-705 Instrument lubricant, 1 gallon

Use self-seal pouches to autoclave and store your instruments.

#164-SS1 Self-seal clear pouch, 3.5" x 8", box of 200

#164-SS2 Self-seal clear pouch, 3.5" x 10", box of 200

#164-SS4 Self-seal clear pouch, 5" x 10.5", box of 200

(Other sizes available)



Use chemical or biological indicators regularly to ensure your autoclave is working properly.

#800-1243B Steam sterilization integrator, bag of 100

#800-1262P Attest® biological indicators, box of 25

#800-116 Attest® biological indicator incubator

For instruments that can't be autoclaved or only require High Level Disinfection (HLD), use a chemosterilant soak.

Check with your Territory Manager to see what cold soak is right for you.



### Aldehydes

#009-28+ BM-28 PLUS™ activated 2% glutaraldehyde disinfectant, 4L

#517-11-2801 MetriCide™ 28 high-level disinfectant, 1 gallon

#421-20394 Cidex OPA, 1 gallon  
HLD 5 minutes at 25°C or 12 minutes at 20°C



## Environmental Infection Control

Wipe down all surfaces regularly.

#907-100906722 Accel® PREvention™ 10" x 10" wipes tub of 60 wipes

#517-11-1100 anti-microbial canister, 160 wipes

#507-PO10-IP70 Isopropyl alcohol, clear, 500mL

#507-PO10-IP70 Isopropyl alcohol, clear, 500mL

Don't forget to clean your stethoscope and blood pressure cuffs.



# Effective exudate control.

## Biatain Fiber with HexaLock Technology

The unique HexaLock Technology integrates gelling fibers in a reinforced hexagon net, ensuring minimal shrinkage and minimising the risk of gap creation and exudate pooling.

**Locks in exudate** to reduce exudate pooling and thereby minimise the risk of infection and maceration

**Minimal shrinkage** to minimise risk of gap creation and exudate pooling

**Cohesive gel** for easy removal and minimal risk of leaving residues behind

To learn more, please visit [biatainfiber.com](http://biatainfiber.com)



Learn more about Biatain Fiber  
Scan the QR code



**Biatain® Fiber**







### Sween® Unscented Lotion

Non-occlusive lotion helps keep skin healthy. Excellent massage lubricant. CHG, latex, and nitrile glove compatible.

#310-7222 2oz/60mL bottle  
#310-7223 4oz/120mL bottle  
#310-7224 21oz/620mL bottle



### Sween® 24 Moisturizing Skin Protectant

Once-a-day application that contains 6% dimethicone. Fragrance-, lanolin-, and alcohol-free. CHG compatible.

#310-7094 3oz/90g tube



stevens.ca/qr-msec20-12

Coloplast Skin Care Resources



### Gentle Rain® Extra Mild Shampoo/Skin Cleanser

Contains mild ingredients and extra emollients to gently clean and moisturize fragile, sensitive skin.

#310-7233 600mL  
#310-7234 1000mL



### Sween® Cream

Contains natural vitamins A and D. Soothes red, sore, dry skin. Absorbs quickly into the skin and hence does not interfere with the adhesion of tape.

#310-7067 3oz/90mL tube  
#310-7068 6.5oz/190mL tube



### Atrac-Tain® Moisturizing Cream

10% Urea and 4% AHA softens and exfoliates severely dry, scaly, calloused, cracked skin. Fragrance- and preservative-free. Latex, nitrile and vinyl glove compatible.

#310-1814 50z/140mL tube



#800-CS102 #800-CS103



### Adhesive Dressings

#002-P150115 1 3/4" x 2", fingertip  
#001-SC7619P 3/4" x 3" strips  
#800-CS102 3/4" x 3" strips  
#800-CS103 1 7/8" x 4", knee/elbow  
#001-SCSPOT 1" dia. spot  
#195-7646013 5cm x 7.2cm fabric patch

box of 100  
box of 100  
box of 80  
box of 10  
box of 100  
box of 100



### Leukoplast® Elastic Adhesive Bandages

#195-7645617 1.9cm x 7.6cm, strips  
#195-7645606 2.5cm x 7.6cm, strips  
#195-7645616 2.2cm, round spots  
#195-7645613 3.1cm x 2.5cm, oval spots  
#195-7645612 3.8cm x 5cm, regular patches  
#195-7645610 5cm x 7.6cm, large patches  
#195-7645607 3.8cm x 7.6cm, knuckle  
#195-7645618 6.3cm x 3.5cm, toe shield (no pad)  
#195-7645608 5.4cm x 3.8cm, small fingertip  
#195-7645609 5cm x 6.3cm, large fingertip  
#195-7645611 7.6cm x 7.6cm, four-wing

box of 100  
box of 100  
box of 100  
box of 100  
box of 100  
box of 50  
box of 100  
box of 100  
box of 100  
box of 50  
box of 50



### Conforming Bandages

#355-A205-X 2" x 4.1yd box of 12  
#355-A305-X 3" x 4.1yd box of 12  
#355-A405-X 4" x 4.1yd box of 12  
#355-A605-X 6" x 4.1yd box of 6



### Elastic Bandages

#379-Q233/50 2" x 5 yards  
#379-Q233/75 3" x 5 yards  
#379-Q233/100 4" x 5 yards  
#379-Q233/150 6" x 5 yards

bag of 12

#379-Q233/75 #379-Q233/100



#379-Q233/150 #379-Q233/50



### Tegaderm™ Transparent Film Dressing

Transparent film dressings designed for protecting skin and wound sites. Frame style allows customization of shape and size to fit any site. Unique sacral shape is designed for challenging and difficult-to-dress areas.

#800-1624W 6cm x 7cm, frame style box of 100



### Coban™ Self-Adherent Wrap

Cohesive material makes it stick to itself but not to other materials or skin.

#800-1581-X 1" x 5yd  
#800-1582-X 2" x 5yd  
#800-1583-X 3" x 5yd  
#800-2084 4" x 5yd  
#800-1586 6" x 5yd



### Hypafix® Conformable Adhesive Retention Tape

- Self-adhesive, non-woven fabric for dressing retention
- Permeable to both air and moisture, helping to reduce the risk of maceration
- Radiotransparent, conformable and can be cut to required size

#195-7144302 10cm x 10m  
#195-7144301 5cm x 10m





### Non-Woven Gauze Sponges, Non-Sterile

- Increased patient comfort with soft viscose/poly blend non-woven material
- No chemicals added—no allergic reactions
- Economical alternative to cotton gauze sponges
- 4-ply

#355-A2101-CH-X 2" x 2" (5cm x 5cm)  
 #355-A2103-CH-X 3" x 3" (7.5cm x 7.5cm)  
 #355-A2100-CH-X 4" x 4" (10cm x 10cm)

**sleeve of 200**  
**sleeve of 200**  
**sleeve of 200**



### All-Gauze Sponges, Non-Sterile

- 100% pure white cotton
- No chemical additives
- For multipurpose applications

#355-B3001-X 8-ply, 2" x 2" (5cm x 5cm)  
 #355-B3003-X 8-ply, 3" x 3" (7.5cm x 7.5cm)  
 #355-B3005-X 8-ply, 4" x 4" (10cm x 10cm)  
 #355-B3002-X 12-ply, 2" x 2" (5cm x 5cm)  
 #355-B3004-X 12-ply, 3" x 3" (7.5cm x 7.5cm)  
 #355-B3006-X 12-ply, 4" x 4" (10cm x 10cm)  
 #355-B3009-X 12-ply, 8" x 4" (20cm x 10cm)

**sleeve of 200**



### Curity™ AMD Antimicrobial Packing Strips

Impregnated with PHMB (Polyhexamethylene Biguanide), a powerful, yet safe antiseptic that has a broad-range effectiveness against gram positive and gram negative microorganisms including some multi-drug resistant strains such as MRSA.

#102-7831AMD-X ¼" x 1yd  
 #102-7832AMD-X ½" x 1yd  
 #102-7833AMD-X 1" x 1yd



### Steri-Strip™ Adhesive Skin Closures, Sterile (Reinforced)

Strips of porous, non-woven backing material, coated with a pressure-sensitive, hypoallergenic adhesive and reinforced with filaments for added strength.

#800-R1540 ⅛" x 3", 5 strips/envelope  
 #800-R1542 ¼" x 1½", 6 strips/envelope  
 #800-R1541 ¼" x 3", 3 strips/envelope  
 #800-R1546 ¼" x 4", 10 strips/envelope  
 #800-R1547 ½" x 4", 6 strips/envelope

**box of 50 envelopes**



### Dukal Premium Non-Woven Sponges

One of the thickest and softest on the market. The heavier fabric offers quick wicking for maximum performance, which results in a reduction of the number of sponges used. They are ideal for delicate procedures and sensitive skin.

#745-6112 4-ply, 2" x 2", non-sterile  
 #745-6113 4-ply, 3" x 3", non-sterile  
 #745-6114 4-ply, 4" x 4", non-sterile

**pack of 200**



### Non-Woven Sponges, 2's

Increased comfort with highly absorbent, soft, non-linting, breathable non-woven material. Two sponges per envelope.

#355-A2240-X 2" x 2", 4-ply  
 #355-A3340-X 3" x 3", 4-ply  
 #355-A4440-X 4" x 4", 4-ply

**tray of 50 envelopes**



### General Purpose Non-Woven Sponges

Made of a soft viscose/poly blend, these sponges are highly absorbent, non-woven, and non-linting.

#680-2100-CH 4-ply, 4" x 4"  
**pack of 200**



### Gauze Sponges, 1's

- Highly absorbent
- Provides optimal protection and comfort
- Heat-sealed envelope to insure safety

#355-B2212-X 2" x 2", 12-ply  
 #355-B3312-X 3" x 3", 12-ply  
 #355-B4412-X 4" x 4", 12-ply

**tray of 100**



### Kerlix™ AMD Antimicrobial Gauze Bandage Rolls

With antimicrobial capabilities, Kerlix™ AMD rolls provide fast-wicking action, superior aeration, and excellent absorbency.

#102-3332 6-ply, 4.5" x 4.1yd



### Non-Adherent Pads

- Absorbent
- Protects without sticking to the wound

#002-P150155 2" x 3"  
 #002-P150160 3" x 4"

**box of 100**



### Dukal™ Wound Closure Strips

- Excellent alternative for sutures or staples
- Each strip consists of thin, reinforced non-woven fibres that provide incision support to maintain closure
- Translucent and breathable

#745-5158 1" x 4", 4 strips per pack  
 #745-5157 ½" x 4", 6 strips per pack  
 #745-5152 ¼" x 1½", 6 strips per pack  
 #745-5151 ¼" x 3", 3 strips per pack  
 #745-5156 ¼" x 4", 10 strips per pack  
 #745-5150 ⅛" x 3", 5 strips per pack

**box of 50 packs**





### Mepore® Dressing

An air-permeable, self-adhesive and absorbent dressing designed especially for surgical wounds. The skin-friendly, water-based, solvent-free polyacrylate adhesive provides a secure yet gentle fixation.

#396-670800 6cm x 7cm **box of 60**  
 #396-670900 10cm x 9cm **box of 50**  
 #396-671000 15cm x 9cm **box of 50**  
 #396-671200 25cm x 9cm **box of 30**



## Adhesive Surgical Tapes

**3M**

### Blenderm™ Plastic Tape, Occlusive, Transparent

#800-1525-1 1" x 5yd **box of 12 rolls**  
 #800-1525-2 2" x 5yd **box of 6 rolls**



**3M**

### Durapore™ Ribbon Tape

#800-1538-1/2 1/2" x 10yd **box of 24 rolls**  
 #800-1538-1 1" x 10yd **box of 12 rolls**  
 #800-1538-2 2" x 10yd **box of 6 rolls**  
 #800-1538-3 3" x 10yd **box of 4 rolls**



**3M**

### Transpore™ Porous Tape, Transparent

#800-1527-1/2 1/2" x 10yd **box of 24 rolls**  
 #800-1527-1 1" x 10yd **box of 12 rolls**  
 #800-1527-2 2" x 10yd **box of 6 rolls**  
 #800-1527-3 3" x 10yd **box of 4 rolls**



**3M**

### Transpore™ Tape, White

#800-1534-1/2 1/2" x 10yd **box of 24 rolls**  
 #800-1534-1 1" x 10yd **box of 12 rolls**  
 #800-1534-2 2" x 10yd **box of 6 rolls**



**3M**

### Micropore™ Paper Tape, White

#800-1530-1/2 1/2" x 10yd **box of 24**  
 #800-1530-1 1" x 10yd **box of 12**  
 #800-1530-2 2" x 10yd **box of 6**  
 #800-1530-3 3" x 10yd **box of 4**



**3M**

### Kind Removal Silicone Tape

#800-2770S-1 1" x 54" **box of 100 rolls**  
 #800-2770S-2 2" x 54" **box of 50 rolls**  
 #800-2770-1 1" x 5.5yd **box of 12 rolls**  
 #800-2770-2 2" x 5.5yd **box of 6 rolls**



### Stanhexidine® 2% w/v Aqueous Chlorhexidine Gluconate Solution

- Contains 4% v/v isopropyl alcohol
- Anti-bacterial
- 450mL flip-cap bottle

#432-L0000009



### Alcohol Prep Pads

70% v/v isopropyl alcohol.

#001-SC-APL large **box of 100**

#001-SC-APM medium **box of 200**



### Benzalkonium Chloride Antiseptic Towelettes

- 200mm X 127mm
- 0.13% w/w Benzalkonium Chloride (applied to the skin)

#821-126-01 **box of 100**



### Isopropyl Alcohol

- 70%

#507-PO10-IP70 500ml bottle clear  
 #680-30070 4L jug



### Tech-Med® Alcohol Bottle

- Unique push down dispensing design prevents cross contamination by not allowing dispensed fluid to return to bottle
- Locking cap prevents fluid from leaking out
- 9oz, imprinted

#745-4024



### Betadine® Solution

10% povidone-iodine topical solution USP. 500mL

- Film forming
- Non-staining
- Water soluble

#036-08617



### Hydrogen Peroxide

- 3% v/v topical solution USP
- 500mL bottle

#248-SHP-500



### Large Solu-Iv Swabstick

2% CHG w/ 70% Isopropyl Alcohol

#800-S10203-BX **box of 50**







### Lactated Ringer's Injections, USP

Comes in a non-PVC, non-DEHP Excel® container. Constructed from a multilayered film specifically developed for parenteral therapies, the container's solution contact layer is a rubberized copolymer of ethylene and propylene that contains no plasticizers and exhibits virtually no leachability.

#971-L7502-X 250mL bag  
#971-L7501 500mL bag  
#971-L7500 1000mL bag



### 0.9% Sodium Chloride Injections, USP

- Available in a PAB® bag or Excel® bag
- Facilitates inspection of admixtures
- Injections supplied in a PAB® bag do not require any overwrap

#971-L8002 250mL Excel® bag  
#971-L8001 500mL Excel® bag  
#971-L8000 1000mL Excel® bag

#971-S8004-526400-X  
100mL fill volume in a 150mL PAB® bag



### Saline Solution

0.9% sodium chloride (NaCl) solution in partial additive bag.

#971-S8004-526400  
100mL/150mL **case of 64**



### 0.9% Sodium Chloride Irrigation

500mL irrigation bottle, screw cap.  
#372-6280



### 3.3% Dextrose and 0.3% Saline IV Solution

#822-JB1034 1000mL bag  
**case of 12**



### 5% Dextrose IV Solution

#971-L5102 250mL bag  
#971-L5101 500mL bag  
#971-L5100 1000mL bag



### B. Braun® Introcan Safety Catheter

#971-4252519-02 22G x 1"  
**case of 200**



### Surflo® Winged Infusion Set

23G x 3/4", 12" tubing, ultra thin wall.  
#691-SV23BLK **box of 50**



### Sterile Water

#971-L8502 250mL bag for IV  
#971-L8501-01 500mL bag for IV  
#971-L8500 1000mL bag for IV  
#010-07990250 500mL injectable  
#372-6290 500mL bottle  
#372-6291 1000mL bottle



### Insyte® Autoguard™ Shielded IV Catheter

- Made of BD Vialon™ biomaterial
- Reduces IV restarts due to reduced complications
- With the push of a button, the needle instantly retracts reducing the risk of accidental needlestick injuries

#333-381423 22G x 1", 35mL/min  
#333-381412 24G x 0.75", 20mL/min  
**box of 50**



### Saf-T-Intima™ Closed IV Catheter System

- Made of BD Vialon biomaterial
- Has wings and extension tubing
- Removable PRN adapter

#333-383318 22mL/min, 24G x 3/4"  
**box of 50**



### Angiocath™ Peripheral Venous Catheter

- Catheter gauge: 24 G, 0.7 mm
- Catheter length: 3/4", 19 mm
- Flow rate: 17 mL/min
- Material: FEP polymer

#333-381112 **box of 50**



### Insyte® Autoguard™ Non-Winged BC Shielded IV Catheter

- Blood-control technology proven to reduce the risk of blood exposure by 95%
- BD Instaflash™ needle technology confirms vessel entry to reduce hit-and-miss insertion providing greater first-stick success for difficult or compromised veins

#333-382512 24G x 3/4" **box of 50**  
#333-382523 22G x 1" **box of 50**



### MicroClave® Clear Neutral Displacement Connector

- Effectively clears blood and blood residual with low flush volumes
- Does not require a change in clinical practice or technique

#179-MC100 **case of 100**







### SOL-M™ Blood Collection Tube Holder

#348-110201050001

case of 4 boxes of 250



### SOL-M™ Tube Holder with Luer Lock Adapter

#348-110201050003

case of 4 boxes of 200



### SOL-CARE™ Safety Winged Infusion Set

DEHP-free tubing is soft, flexible and kink-resistant

#348-213412SI 21G x 3/4" x 12"

#348-233412SI 23G x 3/4" x 12"

#348-253412SI 25G x 3/4" x 12"

case of 4 boxes of 50



### SOL-CARE™ Safety Blood Collection Needle with Luer Adapter

A safety-engineered device intended for venous blood collection. The safety mechanism can be activated immediately after blood draw, helping to protect you against needlestick injury.

#348-110201030001 21G x 3/4" x 12" box of 50

#348-110201030002 21G x 3/4" x 12", w/holder box of 25

#348-110201030005 23G x 3/4" x 12" box of 50

#348-110201030006 23G x 3/4" x 12", w/holder box of 25

#348-110201030009 25G x 3/4" x 12" box of 50

#348-110201030010 25G x 3/4" x 12", w/holder box of 25



### SOL-CARE™ Safety Multi-Sample Blood Collection Needle

- Safety-engineered tool for multi-sample venous blood collection
- The device is designed with a safety mechanism covering the needle after use
- In the activated position, the needle cover guards against accidental needlestick during normal handling and disposal

#348-SMSN22125 22G x 1 1/4" case of 10 boxes of 48

#348-SMSN21125 21G x 1 1/4" case of 10 boxes of 48

#348-SMSNH22125 22G x 1 1/4", with holder case of 8 boxes of 50

#348-SMSNH21125 21G x 1 1/4", with holder case of 8 boxes of 50



**B. BRAUN** SHARING EXPERTISE **LATEX-FREE**

### Primary Gravity IV Set

- 15 drops/mL
- Universal spike
- Slide clamp
- Backcheck valve
- CARESITE® valve injection sites 87" and 12" above distal end
- Roller clamp
- SPIN-LOCK® connector
- DEHP-free
- Priming volume: 20mL
- Length: 104" (264cm)

#971-354203 **case of 50**



**B. BRAUN** SHARING EXPERTISE **LATEX-FREE**

### Primary IV Administration Set

- 15 drops/mL
- Universal spike
- Roller clamp
- Backcheck valve
- Injection sites 6" and 60" above distal end
- Spin-Lock® connector
- Priming volume: 17mL
- Length: 84" (213.4cm)
- DEHP-free

#971-V1443

**case of 50**



**B. BRAUN** SHARING EXPERTISE **LATEX-FREE**

### Primary IV Administration Set

- 15 drops/mL
- Universal spike
- Roller clamp
- Slide clamp
- 1 injection site 26" above distal end
- Spin-Lock® connector
- Priming volume: 18mL
- Length: 93" (236cm)
- DEHP-free

#971-V1415-15

#971-V1415-15-X

**case of 50 each**



**B. BRAUN** SHARING EXPERTISE **LATEX-FREE**

### Primary IV Administration Set

- 15 drops/mL
- Universal spike
- Roller clamp
- Slide clamp
- No injection sites
- Spin-Lock® connector
- Priming volume: 18mL
- Length: 85" (215.9cm)
- DEHP-free

#971-352601



**B. BRAUN** SHARING EXPERTISE **LATEX-FREE**

### Primary Basic IV Administration Set

- 15 drops/mL
- Universal spike
- Roller clamp
- Injection sites: 28" and 6" above distal end
- Spin-Lock® connector
- Priming volume: 18mL
- Length: 84" (213.4cm)
- DEHP-free

#971-V1421-20

#971-V1421-20-X

**case of 50 each**



**B. BRAUN** SHARING EXPERTISE **LATEX-FREE**

### Filterflow® Filtered Administration Set with 0.2 micron Filter

- 15 drops/mL
- Universal spike
- 0.2 micron filter
- Roller clamp
- Slide clamp
- 1 injection site
- Spin-lock® connector
- Priming volume: 21mL
- Length: 91" (226cm)
- DEHP-free
- Latex-free

#971-352424 **case of 50**



**B. BRAUN** SHARING EXPERTISE **LATEX-FREE**

### IV Administration Set

- 15 drops/mL
- Universal spike
- Roller clamp
- Injection site 6" above distal end
- Spin-Lock® connector
- Priming volume: 16mL
- Length: 81" (205.7cm)
- DEHP-free

#971-V1390

#971-V1390-X **case of 50 each**



**B. BRAUN** SHARING EXPERTISE **LATEX-FREE**

### IV Administration Set

- 15 drops/mL
- Universal spike
- Roller clamp
- Slide clamp
- Spin-Lock® connector
- Priming volume: 14mL
- Length: 79" (201cm)
- DEHP-free

#971-V1402

**case of 50**



**B. BRAUN** SHARING EXPERTISE **LATEX-FREE**

### IV Administration Set

- 15 drops/mL
- Universal spike
- Roller clamp
- Spin-Lock® connector
- "Y"-site 6" above distal end
- Priming volume: 14mL
- Length: 73" (185.4cm)
- DEHP-free

#971-V1416

**case of 50**



# DESIGNED FOR UNEXPECTED HEROES

## Are you prepared for sudden cardiac arrest emergencies?

ZOLL Automatic External Defibrillators (AEDs) empower rescuers by delivering a defibrillating shock and unparalleled technology that provides rescuers with real-time, integrated CPR feedback for both rate and depth of compressions.

### HOPE IS IN YOUR HANDS

To learn more about the ZOLL AED Plus  
contact your Stevens Medical representative.

**ZOLL**



## ZOLL

### AED Plus® Automated External Defibrillator

Features Real CPR Help®: a tool that is able to actually see what you are doing and provide feedback to help you do it well. Audio and visual prompts help you rescue with confidence and clarity, unmatched by any other AED (Automated External Defibrillator).

- Not pushing hard enough? It will tell you when to push harder
- Not pushing fast enough? A metronome will lead you to the right rate
- Zoll AED will even show you the depth of each compression — in real time

#724-22600000102011060

#724-22600000101011260

#724-22600710702011060

#724-22600710701011260

#724-8900-0800-01

#724-8000-0855

semi-automatic AED, English

semi-automatic AED, French

fully-automatic AED, English

fully automatic AED, French

spare pads

cabinet



MCN CP 1705 0039



## 911 Doctor's Emergency Kit

Kit contains:

- 1 x "D" cylinder w/o oxygen
- 1 x Clamshell carrying bag
- 8 x Airways (40mm-110mm in a ZipLoc package)
- 1 x Disposable adult BVM
- 1 x 0-15LPM adjustable regulator
- 1 x Wrench
- 1 x Adult high concentration mask w/ tubing
- 1 x Pediatric high concentration mask w/ tubing
- 1 x Adult cannula
- 1 x Pediatric cannula
- 2 x Pairs of latex-free gloves
- 1 x Surgitube tubular gauze
- 1 x CPR Kit
- 1 x MDI spacer
- 1 x Pocket vent mask with oxygen concentrator

#991-911-X



stevens.ca/qr-msec20-18

Zoll AED Plus resources

## CPR-Uni-Padz

- Universal Adult or Child Electrode Pads
- Chest sensor includes Real CPR Help® real-time feedback for depth and rate of compressions during CPR
- 5-year shelf life
- Rescue kit: includes gloves, scissors, razor, barrier mask, wet/dry wipe
- Compatible with ZOLL AED 3 and ZOLL AED 3 BLS

#724-8900-000260



## Crash Cart Kit

Kit contains:

- 1 x "D" cylinder w/o oxygen
- 1 x 0-15LPM adjustable regulator
- 1 x Wrench
- 1 x Adult high concentration mask w/ tubing
- 1 x Pediatric high concentration mask w/ tubing
- 1 x Adult cannula
- 1 x Pediatric cannula
- 1 x Adult BVM
- 1 x Set of airways 40mm-110mm

#991-411-X



## CPR Protective Barrier

With valve, high efficiency filter (99.2%) and exhalation port. Comes in 11cm x 11cm, black key ring pouch. Includes one pair of latex-free gloves and one alcohol antiseptic swab.

#991-3052B-X



## Guedel Oral Airway

90mm, yellow, adult medium

#991-2109

\*Other sizes available



## Emergency Laryngoscope Set

XP with F/O 2.5V handle. Includes blades and hard case.

#799-F-257.10.815



## Guedel Oral Airways

- Latex-free, disposable, non-sterile
- Colour coded
- Individually wrapped in an easy-open packaging
- Set of 8 airways in a plastic case (sizes: 40mm, 50mm, 60mm, 70mm, 80mm, 90mm, 100mm, and 110mm)

#991-2200-10 kit

**These airways are also available for sale individually!**



## Disposable Resuscitator Mask

Oxygen reservoir system, O<sub>2</sub> tubing with pop-off.

#991-4010F child w/filter (10-30kg/22-66lb)

#991-4025F infant w/filter (<10kg/22lb)

#991-4050F adult w/filter (30kg/66lb+)

#991-4010 child (10-30kg/22-66lb)

#991-4025 infant (<10kg/22lb)

#991-4050 adult (30kg/66lb+)



## Adult Straight Tip Cannula

- 7ft. (2.1m) SureFlow™ tubing
- Soft curved tip

#001-SC-OT-1012SC



## Aerosol Kit

- Single use

Includes:

- Aerosol mask
- Updraft #1800 nebulizer
- With 7ft (2.1m) SureFlow™ tubing

#991-1807 with pediatric mask

#001-SC-OT-1805 with adult mask



## Oxygen Kit

Elongated, medium-concentration adult mask with 7ft (2.1m) SureFlow™ tubing.

#001-SC-OT-1041



## MDI Rigid Spacer

180cc

#991-1400S with pediatric mask

#991-1400M with child/small adult mask

#991-1400L with adult mask



## Suction Connecting Tubing

#184-70-8072 1/4" x 72"

#184-70-8120 1/4" x 120"

**case of 50**

**case of 25**



## CPR Barrier

CPR training print face shield with high-efficiency filter (99.2%).

#991-3000

**case of 100**



## CPR Pocket Ventilator

Includes head strap, oxygen inlet, filter, and clamshell hard case.

#991-3100-50







**Kimberly-Clark**  
Kleenex® Ultra Folded  
Paper Towels

- 2-ply
- White
- Size: 16.3" x 8.5"  
(41.4cm x 21.6cm)

#151-28791

**case of 30 94-sheet packs**



**Kimberly-Clark**  
Scott® Paper Towels

- M-folded
- 7.8" x 12.4" (19.8cm x 31.5cm)
- 1-ply
- White

#151-01960

**case of 25 175-sheet packs**



**Kimberly-Clark**  
Scott® Paper Towels

- Single-fold
- 10.5" x 9.3" (26.7cm x 23.6cm)
- 1-ply
- White

#151-01700

**case of 16 250-sheet packs**



**Facial Tissue**

2-ply.

#001-SC-FT-1114-X 19cm x 12.7cm

**box of 80 sheets**

#001-SC-FT-1222-X 19cm x 21cm

**box of 100 sheets**



Kleenex® Facial Tissue

White, 2-ply sheets in a pop-up box.

#151-21400 8.4" x 8" sheets

**case of 36 100-sheet boxes**

#151-21400-X 8.4" x 8" sheets

**100-sheet box**



White Swan® Quick Wipes

White, 2-ply, 8.3" x 4.5"  
sheets in a pop-up box.

#340-17900298

**80-sheet box**



**Kimberly-Clark**  
Scott® Essential Plus  
Hard Roll Paper Towels

- 8" x 600' (20.3cm x 182.9m) rolls
- White

#151-50606

**case of 6 rolls**



**Paper Towels**

- 3-ply
- White, embossed tissue
- Size: 13½" x 18"

#001-SC054

**case of 500**



**Kimberly-Clark**  
Scott® Bathroom Tissue

#151-48040 2-ply, white,  
4.1" x 4.0" sheets, standard rolls

**case of 40 550-sheet rolls**

#151-07805 2-ply, white,  
3.55" x 1000ft jumbo rolls, non-wrapped

**case of 12 rolls**



**Kimberly-Clark**  
Scott® Paper Towels

- Multi-fold
- 9.5" x 9.25" (24cm x 23.5cm)
- 1-ply, white

#151-01804

**case of 16 250-sheet packs**



**Kimberly-Clark**  
Windows® Universal Folded  
Towel Dispenser

- White
- Size: 13.31" x 18.85" x 5.85"  
(33.8cm x 47.9cm x 14.9cm)

#151-09906



**Kimberly-Clark**  
Scott® Kitchen Roll Towels

- Perforated rolls with  
absorbency pockets
- 11" x 8.78" / 27.9cm x 22.3cm
- 1-ply, white

#151-41482

**case of 20 128-sheet rolls**

## Tork Elevation Mini Xpress™ Hand Towel Interfold Dispenser

- Measuring just 12" high, the Mini Xpress™ is perfect for wall areas where space is at a premium
- The one-at-a-time dispensing means you touch only the towel you use—no cross-contamination from levers or cranks. Plus the hygiene shield protects your towel until used
- Easy and intuitive maintenance, refilling, and use

#340-141014241



Looking for a dispenser?



Your local Stevens Rep can help!



### Prime Source® Single-Fold Towels

- High-quality, value-priced paper towels
- Ideal for all public environments
- Size: 10.5" x 9.3"

#340-75000250 natural  
#340-75000251 white

**case of 16 x 250**



### Multi-Fold Towels

- High-quality, value-priced paper towels
- Ideal for all public environments

#340-14100538 white, universal

**case of 16 x 250**

#340-75000254 kraft

**case of 16 x 250**

### Tork Premium Hand Towel Interfold

- The unique 4 panel interfolded system creates ready-to-use towels, with no towel "fall out"
- Portion control delivers one towel every time
- Easy to refill

#340-14100979 **case of 16 x 94**



### Tork Perforated Roll Towels

- High-quality, value-priced paper towels
- White
- 84 sheets per roll

#340-14100122

**case of 30 rolls**

### Tork Advanced Perforated Towel Roll

- 120 towels per roll
- Towel size: 11" x 6¾"
- Colour: white

#340-14100107

**case of 30 rolls**

### Tork Bathroom Tissue

- White, 2-ply
- 10.1cm x 9.5cm sheets
- 500 sheets per roll

#340-14106131

**case of 48 rolls**

### Prime Source® Roll Towels

- High-quality, value-priced paper towels
- Ideal for all public environments
- Fits 84TR, 86TR, 87TR dispensers

#340-75000257 8" x 350ft, natural

#340-75000256 8" x 350ft, white

**case of 12 rolls**



### Examination Table Paper

- #001-SC356 smooth, 18" x 225ft rolls
- #001-SC355 crepe, 18" x 125ft rolls
- #001-SC918 smooth, 18" x 225ft rolls
- #001-SC818 crepe, 18" x 125ft rolls
- #001-20210 smooth, 17¾" x 210ft rolls
- #001-SC921 smooth, 21" X 225ft rolls
- #001-00110 crepe, 17¾" x 135ft rolls
- #884-02775 lite crepe, 17¾" x 248ft rolls

#### Western Canada only:

- #001-00201 smooth, 17" x 248ft rolls
- #001-00210 smooth, 17¾" x 248ft rolls
- #001-00101 crepe, 17" x 135ft rolls
- #884-30101 premium crepe, 17" x 248ft rolls

**case of 12 rolls**



### MediShorts® Disposable Exam Shorts

Made of durable non-woven material with latex-free waist.

- #200-10000 18"-44", small / medium
- #200-10001 22"-48", large / X-large

**case of 50**



### Disposable Pillowcase

- Tissue/Poly construction protects pillow from moisture and cross contamination
- Size: 21" x 30"; fits standard size pillows
- White

#002-P230023

**case of 100**



### Disposable Pillowcase

White, 21" x 30"

#200-360

**case of 100**



### Vinyl Pillowcase

- Zipper closure
- Fits standard size pillow
- Latex-free
- PVC. 21" x 27".
- Thickness - 4MIL (.10MM)

#139-3857

**each**



### Disposable Exam Gowns

- Size: 30" x 42"
- Tissue/poly/tissue
- Blue

#001-SC731

**case of 50**



### Disposable Exam Capes

- Size: 30" x 21"
- Tissue/poly/tissue
- Blue

#002-P750003

**case of 100**



### Disposable Exam Gown

- Size: 30" x 42"
- Knee-length coverage
- 3-ply
- White

#002-P750333

**case of 50**



### Disposable Exam Gowns

- Size: 28" x 42"
- Layers: tissue/poly/tissue
- Blue

#884-22842

**case of 100**



### Isolation Gown

- Yellow with synthetic elastic cuff.
- #355-D-8009 regular

#680-8010

X-large

**case of 50**



### Impervious Isolation Gown

- Blue, regular with thumb loop.

#355-8025

**case of 75**

## Disposable Drapes

- 2-ply tissue
- White
- Size: 40" x 48"

#001-SC400

**case of 100**



✓ LATEX-FREE

## Caliber™ Surgical Drapes, 3-ply

- Excellent for use during minor surgical procedures
- Two layers of absorbent tissue with a polyethylene lining in between help eliminate risk of bacterial migration and moisture strike through
- Size: 18" x 26"

#745-20-002

non-fenestrated

#745-20-001

2.75" fenestration

**box of 50**



## Disposable Drapes

- Fanfold, 2-ply tissue
- White
- Size: 36" x 48"

#001-SC444

**case of 100**



## Disposable Drapes, White

- Fanfold
- 2-ply tissue
- Size: 36" x 40"

#200-50363

**case of 100**



## Emesis Basin

- Capacity: 24oz
- Graduated: 100cc increments from 100cc to 840cc, as well as every 4oz from 4oz to 24oz
- Sterilizable
- Made of plastic
- Colour: blue

#193-00063



## Attends

### Dri-Sorb® Underpads

- Light absorbency level
- Perfect for standard, everyday use
- Size: 17" x 24"

#850-31925

**case of 30 bags of 10**



## Prevail® Underpads

- Designed to protect beds, chairs, and other surfaces
- Fluff, green, 23" x 36"

#740-UP-150

**case of 10 packs of 15**



## Paper Food Trays

- Natural pulp fibre
- Beige

#584-071126-X

7 $\frac{3}{8}$ " x 5 $\frac{1}{2}$ " x 1 $\frac{1}{2}$ "

#584-071133-X

8.4" x 6.3" x 1.6"



## Waxed Paper Cups

#001-SC-PC-04 4oz

#001-SC-PC-05 5oz

**pack of 100**



## Medicine Cups

- Clear poly
- 30mL capacity
- Graduated in tbs, fl. oz., drams, cc, and mL

#001-SC-MC-01-X

**pack of 100**



## Styrofoam Cups

- Maintain beverages at their optimal temperatures longer
- 7oz
- White

#584-064002

**case of 1000**



## Pro-Medix Paper Portion Cups

- 1oz
- Single use

#001-SC-6400-X

**box of 250**

#001-SC-6400

**case of 5000**



## Pro-Medix Plastic Cups

- Single use
- #001-SC-6210  
5oz, rose

**case of 1000**

#001-SC-6205

5oz, blue

**case of 1000**







### Set of Sundry Jars

- Clear flint glass with overlapping, stainless steel cover
- 7"H x 4 1/4" diameter
- 139-3454: permanently imprinted in blue: "Bandages", "Applicators", "Tongue Depressors", "Gauze", and "Cotton"

#139-3454 labeled **set of 5**  
 #139-3458 unlabeled **each**

### Garbage Bags

#455-2952-10 regular, 20" x 22", white  
 #455-2752-01 regular, 20" x 22", black  
**case of 500**

#455-2771-01 strong, 26" x 36", black  
 #455-2773-01 strong, 30" x 38", black  
**case of 200**

#455-2775-01 strong, 35" x 50", black  
**case of 125**



### Odorgon® Odour Eliminator

- Non-toxic
- Biodegradable
- Environmentally friendly

#741-988-0144 Original, 114mL bottle  
 #741-988-0074 Original, 237mL bottle



### Distilled Water

4L jug.

#320-72404  
 Western Canada, use code: #781-AQ00334  
**case of 4**



### Cotton Balls

Economy, medium, 0.22g

#589-21176

**bag of 2000**



### PSI FLEX-GEL Reusable Cold Packs

- Non-caustic and non-toxic
- Reduces bleeding, hemorrhaging and swelling
- Controls muscle spasms and pain - promotes healing

#061-3200



✓ LATEX-FREE

### Rapid Relief® Cold/Hot Gel Compress with Form-Fit Gel

- Reusable
- Microwaveable
- Non-toxic blue gel

#161-12259-24 5 1/4" x 9"  
 #161-12246-24 4" x 6"

**each**  
**case of 24**



✓ LATEX-FREE

### Plastic Teaspoons

- Disposable
- Single use

#340-75002494

**case of 1000**



✓ LATEX-FREE

### Pro-Medix Straws

- Single use

#001-SC-7000-X

unwrapped

**box of 500**

#001-SC-7000

unwrapped

**case of 20 x 500**



### MedPro Flexible Straws

MedPro Flexible Straws, Wrapped 19 cm x 6 mm  
 500 per box, 20 boxes per case

#635-018-550 wrapped

**500 / box, 20 boxes / case**



Medicom

### Triangular Bandages

- Size: 40" x 40" x 56"
- Cotton-muslin blend
- 2 safety pins included

#680-5460



### Fluid Disposal Kit

Nitrile exam gloves; 124.5cm disposable apron, biohazard pickup scoop/scrapper, CaviWipes disinfectant, surgical face mask w/shield + tie strings, identification tag for universal biohazard kit, infectious waste bag, two Paws antimicrobial hand towelettes, 21g Red Z package, wiping towel, and white twist tie  
 #390-607923ZC



✓ LATEX-FREE

### Rapid Relief® Instant Cold Packs

- Easy activation with new FTP™ technology
- Convenience and portability with no more worries about leaky, messy ice bags
- Fast relief for: bruises and swelling, muscle spasms and pain, headaches and minor injuries

#161-31246-50 4" x 6"

#161-31259-24 5" x 9"



INNOVATEK

✓ LATEX-FREE

### Cotton-Tipped Applicators

- Wooden base

#553-4202-40-100-X 3"

#001-SC-CTA-6-X 6"

**bag of 100**



AMD

### OB-GYN/Throat Swabs

- Rayon tip
- 8" paper shaft

#355-57600


**box of 50**





# STEVENS

*Inspired by the care you deliver*



## Contact Information

### British Columbia



19051-27<sup>th</sup> Avenue  
Surrey, BC V3Z 5T1

Tel: 604-634-3088

Fax: 604-585-0193

Toll-free: 1-800-565-8444

Email: BCCS@stevens.ca

### Manitoba



475 De Baets Street—Unit 2  
Winnipeg, MB R2J 4M3

Tel: 204-885-9440

Toll-free fax: 1-888-640-8088

Toll-free: 1-800-665-0368

Email: MWCS@stevens.ca

### Atlantic Canada



101 Thornhill Drive—Unit 109  
Dartmouth, NS B3B 1S3

Tel: 902-468-7582

Fax: 902-468-7824

Toll-free: 1-800-565-0765

Email: ACCS@stevens.ca

[www.stevens.ca](http://www.stevens.ca)

### Midwestern Canada



2620-61<sup>st</sup> Avenue S.E.  
Calgary, AB T2C 4V2

Tel: 403-640-2858

Fax: 403-640-2976

Toll-free: 1-800-665-0368

Email: MWCS@stevens.ca

### Ontario—Head Office



425 Railside Drive  
Brampton, ON L7A 0N8

Tel: 905-791-8600

Fax: 905-791-6143

Toll-free: 1-800-268-0184

Email: ONCS@stevens.ca

### Quebec



9390 boul des Sciences  
Anjou, QC H1J 3A9

Tél: 514-352-5652

Télé.: 514-352-9214

Sans frais: 1-855-660-7750

Email: QCSAC@stevens.ca

Regular terms and conditions apply. Errors and omissions excepted.  
All products approved for sale by Health Canada at time of printing.

ID# 1550866491

Published May, 2022



M-CAT-MSE-22